Cochlear™ Nucleus® 8 Aqua+ Coils

User Guide

EN

ID

KO

MS

TH

VI

ZH-TW

Contents

EN	Cochlear™ Nucleus® 8 Aqua+ Coils User Guide	1
ID	Panduan bagi pengguna Koil Aqua+ Nucleus® 8 Cochlear™	31
KO	Cochlear™ Nucleus® 8 Aqua+ 코일 사용 안내서	61
MS	Panduan Pengguna Gegendang Aqua+ Cochlear™ Nucleus® 8	91
TH	คู่มือผู้ใช้คอยล์ Cochlear™ Nucleus® 8 Aqua+	121
VI	Hướng dẫn sử dụng Cochlear™ Nucleus® 8 Aqua+ Coil	151
ZH-TW	Cochlear™ (科利耳®) Nucleus® 8 Aqua+ 線圈使用指南	181

Symbols used in this document

Note: Important information or advice.

Tip: Time saving hint.

Caution: Special care to be taken to ensure safety and effectiveness. Could cause damage to equipment.

Warning: Potential safety hazards and serious adverse reactions. Could cause harm to person.

Contents

Symbols used in this document.....	ii
About this guide	3
Cochlear Nucleus 8 Aqua+ Coil	4
Introduction	6
Intended purpose	6
Indications	7
Contraindications.....	9
Intended patient population.....	11
Benefits.....	11
Intended users	11
Use.....	12
Connect the coil.....	12
Wear	13
Wear your sound processor	13
Care	15
Aqua+ coil magnet instructions.....	15
Care and storage.....	16

Troubleshoot.....	17
Warnings.....	19
For parents and carers.....	19
Sound processor and parts	20
Medical treatments	22
Other information	23
Operating Characteristics	23
Product component measurements.....	24
Environmental conditions	24
Summary of safety and clinical performance	24
Labelling symbols	25
Serious incidents.....	27
Electromagnetic compatibility (EMC)	28
Environmental protection	28
Equipment classification.....	28
Legal statement	29
Trademark legal statement.....	29

About this guide

This guide is intended for hearing implant recipients and their carers using the Cochlear™ Nucleus® 8 Aqua+ coils.

This guide can also be used as a reference by hearing care professionals.


Note

Refer to the sound processor user guide relating to the use of the processing unit, batteries and other devices.

Please also refer to your *Important Information* document for essential advice that applies to Cochlear implant systems.

Cochlear Nucleus 8 Aqua+ Coil

1. Cochlear Nucleus 8 Aqua+ Coil
2. Cochlear Magnet
3. Coil connector
4. Coil cable

Device

Cochlear Nucleus 8 Aqua+ Coil

Available in two cable lengths (6 and 8 cm).

Cochlear Nucleus 8 Aqua+ Coil (5(I))

Available in two cable lengths (6 and 8 cm).

Compatible accessories

- Cochlear Aqua+ Compact Sleeve
- Cochlear Aqua+ Power Extend Sleeve
- Cochlear Safety Line
- Cochlear Coil Spacer
- Cochlear Magnet
- Cochlear Magnet Reverse Polarity

Introduction

Intended purpose

Device	Intended purpose
Cochlear Nucleus 8 Aqua+ Coil	The coil is intended to be used in combination with other devices as part of a hearing implant system to supply electrical signals and power between a processing unit and a hearing implant.
Cochlear Nucleus 8 Aqua+ Coil (5(I))	
Cochlear Magnet Cochlear Magnet Reverse Polarity (all strengths, types and colours)	The magnet is intended to align and hold the coil or processing unit to a hearing implant.

Indications

Cochlear Nucleus 8 Aqua+ Coil

The Cochlear Nucleus 8 Aqua+ Coil is indicated for a recipient with a compatible Cochlear Nucleus® implant, a compatible Nucleus processing unit and compatible Cochlear magnet.

Compatible Cochlear Nucleus implants are:

- CI600 Series Implants : CI612, CI622, CI624, CI632
- CI500 Series Implants : CI512, CI513, CI522, CI532, CI551, ABI541
- CI24RE Series Implants: CI422, CI24RE (CA), CI24RE (ST), CI24RE (CS), CI24REH (Hybrid L24), CI8REH (Hybrid S8), CI12REH (Hybrid S12)
- CI24R Series Implants: CI24R (CA), CI24R (CS), CI24R (ST)
- CI24M Series Implants: CI24M, ABI24M, CI 11+11+2M, CI24MH (CI6+16+2M)

Compatible processing unit is the Cochlear Nucleus 8 Processing Unit.

Compatible magnets are Cochlear Magnet strength ½ to 6, ½(I) to 4(I) and Cochlear Magnet Reverse polarity strength ½ to 6.

Cochlear Nucleus 8 Aqua+ Coil (5(I))

The Cochlear Nucleus 8 Aqua+ Coil (5(I)) is indicated for a recipient with a compatible Cochlear Nucleus implant, a compatible Nucleus processing unit and compatible Cochlear magnet.

Compatible Cochlear Nucleus implants are:

- CI600 Series Implants : CI612, CI622, CI624, CI632

Compatible processing unit is the Nucleus 8 Processing Unit.

Compatible magnets are Cochlear Magnet strength 5(I).

Cochlear Magnet

The axially polarised (M) magnet is indicated for use with the Cochlear Nucleus 8 Slimline Coil and Cochlear Nucleus 8 Aqua+ Coil.

The diametric (I) magnet is indicated for use with the Cochlear Nucleus 8 Slimline Coil and Cochlear Nucleus 8 Aqua+ Coil.

Cochlear Magnet Reverse Polarity

The Cochlear Magnet Reverse Polarity is indicated for use with the Cochlear Nucleus 8 Slimline Coil and Cochlear Nucleus 8 Aqua+ Coil.

Cochlear Magnet strength 5(I)

The Cochlear Magnet strength 5(I) is indicated for use with the Cochlear Nucleus 8 Slimline Coil 5(I)) and Cochlear Nucleus 8 Aqua+ Coil (5(I)).

Contraindications

Cochlear Nucleus 8 Aqua+ Coil

The Cochlear Nucleus 8 Aqua+ Coil is not indicated for use with hearing implants or processing units other than those listed in Indications.

The Cochlear Nucleus 8 Aqua+ Coil is not indicated for use with Cochlear Magnet strength 5(I).

Cochlear Nucleus 8 Aqua+ Coil (5(I))

The Cochlear Nucleus 8 Aqua+ Coil (5(I)) is not indicated for use with hearing implants or processing units other than those listed in Indications.

The Cochlear Nucleus 8 Aqua+ Coil (5(I)) is not indicated for use with Cochlear Magnet strength $\frac{1}{2}$ to 6 and $\frac{1}{2}$ (I) to 4(I) and Cochlear Magnet Reverse Polarity strength $\frac{1}{2}$ to 6.

Cochlear Magnet

The Cochlear Magnet is not indicated for use with the Cochlear Nucleus 8 Slimline Coil (5(I)) or Cochlear Nucleus 8 Aqua+ Coil (5(I)).

Cochlear Magnet Reverse Polarity

The Cochlear Magnet Reverse Polarity is not indicated for use with Cochlear Nucleus 8 Slimline Coil (5(I)) or Cochlear Nucleus 8 Aqua+ Coil (5(I)).

Cochlear Magnet strength 5(I)

The Cochlear Magnet strength 5(I) is not indicated for use with the Cochlear Nucleus 8 Slimline Coil and Cochlear Nucleus 8 Aqua+ Coil.

Intended patient population

The sound processor is intended for patients implanted with a compatible Cochlear Nucleus implant. There are no restrictions for the intended patient population of the sound processor and accessories in terms of age, weight, health or other condition.

Benefits

The sound processor operates in conjunction with a compatible Cochlear Nucleus implant system.

Potential benefits of receiving a Cochlear Nucleus implant system include:

- Better understanding of speech in quiet.
- Better understanding of speech in noise.
- Increased satisfaction based on hearing capabilities.

Intended users

Intended users of the Cochlear Nucleus 8 Aqua+ and Cochlear Nucleus 8 Aqua+ (5(I)) coils are:

- Recipients of a compatible Cochlear Nucleus implant.
- Carers of recipients, who carry out recipient functions as needed. Carers may include parents of paediatric recipients who are less than 12 years old, and nurses or other carers of otherwise dependant recipients.
- Audiologists/clinicians.

Use

Connect the coil

The Cochlear Nucleus 8 Aqua+ Coil is used in place of your usual coil while your sound processor is inside your Aqua+. The Aqua+ coil and Aqua+ keep your sound processor dry during use in or around water (e.g. when swimming, surfing, or just having fun splashing around). It has a special plug that creates a seal when connected to your sound processor's coil socket through the coil plug hole in Aqua+. As with your usual coil, a coil magnet screws into the Aqua+ coil and holds it against your implant.

Push the blue Aqua+ coil (with magnet) straight into your sound processor through the hole in Aqua+ until it clicks into place. **Do not twist**.

Note: Refer to your Aqua+ user guide for details about using the Aqua+ including warnings and precautions.

Wear

Wear your sound processor

1. Place the sound processor on your ear, letting the coil dangle.

2. Move the coil sideways and onto your implant.

Indicator lights

Flash of orange every second

What it means

Sound processor flashes while coil is off (or connected to the wrong implant).

People with CI600 Series implants

If you have a CI600 Series implant, avoid sliding the coil sideways onto your implant. This could cause the coil magnet to misalign with your implant. Always place the coil down onto your implant.

To place the coil on your head:


1. Hold the coil slightly above the implant location on your head.
2. Rotate the coil slightly in both directions (clockwise and anti-clockwise).

3. When you feel a strong pull, place the coil on the implant.
4. Rotate the coil to a comfortable position for wearing.

Care

Aqua+ coil magnet instructions


1. Insert the new magnet and turn clockwise until it stops.
2. Turn the magnet a little more until you feel a click.

Tip: The tamper-resistant lock marker aligns with the cable when locked

For Cochlear Magnet strength 5(I) only.

3. Align the cover on the magnet and press until you feel a click.

Care and storage

1. Rinse Aqua+ coil with clean fresh water and dry with a soft cloth.
2. Unplug the blue Aqua+ coil. Grip the black side ribs of the coil connector. Pull the plug straight out.

Caution: Do not pull on the coil cable or twist when you pull it out.

3. Place the blue Aqua+ coil in your dry aid kit overnight.

4. Once they're thoroughly dry, store the Aqua+ coil in the Nucleus 8 Activity Kit's case.

Troubleshoot

Contact your clinician if you have any concerns regarding the operation or safety of your sound processor.

Problem	Resolution
You do not hear sound or sound is intermittent	<ol style="list-style-type: none">1. Make sure the coil cable is fully inserted into the socket on the processing unit.2. Make sure you are using the correct coil magnet for your implant. If unsure, contact your clinician.3. Make sure the coil is oriented properly on your head. See <i>Wear your sound processor</i> on page 134. If you use the Remote Control, turn up the volume.5. If you use the Nucleus Smart App, turn up the volume or sensitivity.6. Try a different program.7. Change the batteries.8. If the problem continues, contact your clinician.

Problem	Resolution
Coil does not attach as strongly as usual	Make sure the coil is oriented properly on your head. See <i>Wear your sound processor</i> on page 13
The sound processor or coil become hot	Remove the processing unit, coil and any cables from your head immediately, disconnect the battery module and contact your clinician.
You experience tightness, discomfort or develop a skin irritation at your implant site	<ol style="list-style-type: none">1. Your coil magnet may be too strong or in contact with your skin. Change to a weaker magnet. See <i>Aqua+ coil magnet instructions</i> on page 152. Try using an adhesive SoftWear pad. For instructions on using SoftWear pads, refer your sound processor user guide.3. If you are using a retention aid, such as a headband, this may be placing pressure on your sound processor. Adjust your retention aid, or try a different aid.4. If the problem continues, contact your clinician.

Warnings

For parents and carers

- Unsupervised use of long cables (e.g. coil or accessory cables) or the headband may present a risk of strangulation.
- Carers must routinely check the device for signs of overheating. Remove the sound processor, coil and any cables immediately if the device becomes hot and contact your clinician.
- Carers must routinely check for signs of discomfort (e.g. if sound is uncomfortably loud). Remove the sound processor immediately if there is any discomfort and contact your clinician.
- Carers must routinely check for signs of discomfort or skin irritation at the implant site. Remove the sound processor immediately if there is any discomfort and contact your clinician.
- Carers must monitor for signs of discomfort or skin irritation if a retention aid (e.g. headband) is used that applies pressure to the sound processor. Remove the aid immediately if there is any discomfort or pain and contact your clinician.
- Sound processors and related accessories contain small parts that alone or in combination may pose a hazard of inhalation, choking or ingestion. Swallowing or inhaling small parts can cause severe or fatal injuries. Use the tamper-resistant locks to keep small parts attached to the processing unit where available, including parts that fit between the processing unit and battery module. Always supervise children under 3 years of age and others who may be at risk of inhalation, choking or ingestion of small parts when they use the sound processor and related accessories. When not in use, keep small parts, and combinations of small parts, out of reach from children. If small parts are swallowed or inhaled, seek immediate medical attention.

Sound processor and parts

- Your sound processor and other parts of the system contain complex electronic parts. These parts are durable but must be treated with care.
- No modification of this equipment is allowed. Warranty will be void if modified.
- Do not place the device or accessories inside any part of your body (e.g. nose, mouth).
- If the coil magnet is too strong or is in contact with the skin, pressure sores may develop at the implant site. If this happens, or if you experience tightness or pain in this area, stop using your sound processor and contact your clinician.
- Do not apply continued pressure to the coil when in contact with the skin (e.g. sleeping while lying on coil, or using tight fitting headwear).
- Turn your sound processor off before cleaning or performing maintenance.
- Do not use cleaning agents or alcohol to clean your sound processor or accessories.

- Do not use your sound processor in an explosive or oxygen-rich environment.
- Do not let your sound processor or parts entangle with any jewellery (e.g. earhook and earrings) or machinery.
- Do not place the sound processor or parts in any household devices (e.g. microwave oven, dryer).
- Do not expose the sound processor or parts to heat (e.g. never leave them in sunlight, behind a window or in a car).
- The magnetic attachment of your sound processor to your implant may be affected by other magnetic sources.
- The magnetic attachment of your sound processor to your implant may affect hearing aids.
- Your sound processor coil and magnet may be affected by metallic or magnetic objects. Keep metallic or magnetic objects away from your coil.
- Store spare coil magnets safely and away from cards that may have a magnetic strip (e.g. credit cards, bus tickets).
- Before activities that create electrostatic discharge (e.g. playing on plastic playground equipment), remove your sound processor. In rare cases, discharge of static electricity can damage or cause your sound processor to shut down. If your sound processor shuts down, it should resume normal operation after restarting it. If static electricity is present (e.g. when putting on clothes over your head, or getting out of a car), before the Cochlear implant system touches any object or person, you should touch something conductive such as a metal door handle.

- Your device contains magnets that should be kept away from life supporting devices (e.g. cardiac pacemakers and ICDs (implantable cardioverter defibrillators) and magnetic ventricular shunts), as the magnets may affect the function of these devices. Keep your sound processor at least 15 cm (6 in) from such devices. Contact the manufacturer of the specific device to find out more.

Medical treatments

Magnetic resonance imaging (MRI)

The processing unit, the Cochlear Nucleus 8 Aqua+ coil, remote control and related accessories (such as the Wireless Programming Pod) are MR Unsafe.

Full MRI safety information is available at www.cochlear.com/warnings or by calling your regional Cochlear office (contact numbers available at the end of this document).

Other information

The coil acts as a transformer coupling that transfers energy and data to the implant.

Materials

Coil	Polypropylene (PP), thermoplastic elastomer (TPE)
Coil magnet casing	Acrylonitrile butadiene styrene (ABS)
Coil cable sheath	Polyvinyl chloride (PVC)
Coil cable plugs	PP and TPE

Operating Characteristics

Coil

Characteristic	Value / Range
Operating voltage	2.0 V to 2.6 V
Operating frequency	5 MHz

Product component measurements

Typical values. The weight is measured with a medium earhook and cable. The weight does not include coil magnet.

Item	Weight	Depth	Diameter
Cochlear Nucleus 8 Aqua+ Coil	3.9 g	6.4 mm	30.7 mm

Environmental conditions

Condition	Minimum	Maximum
Storage and transport temperature	-10° C (+14° F)	+55° C (+131° F)
Storage and transport humidity	0% RH	90% RH
Operating temperature (sound processor)	+5° C (+41° F)	+40° C (+104° F)
Operating relative humidity	0% RH	90% RH
Operating pressure	700 hPa	1060 hPa

The temperature of the sound processor and its accessories may rise by up to 2.7° C (4.86° F) during normal operation, and could result in these components reaching a temperature of +42.7° C (+108.86° F) when operated at the maximum environmental temperature of +40.0° C (+104° F).

Summary of safety and clinical performance

A summary of the safety and clinical performance of the Cochlear Nucleus 8 Aqua+ coils can be found at

<https://ec.europa.eu/tools/eudamed>.

Labelling symbols

The following symbols may appear on your sound processor or accessories and/or packaging:

Refer to instruction manual

Consult instructions for use

Specific warnings or precautions associated with the device, which are not otherwise found on the label

Manufacturer

Compatible sound processors

Compatible implants

Authorised representative in the European Community

Unique Device Identification

Medical Device

Catalogue number

Other information

LOTBatch code

Date of manufacture

Temperature limits

CE₀₁₂₃CE registration mark with notified body number

Rx OnlyBy prescription

Recyclable material

Dispose of electrical components in accordance with your local regulations

Type B applied part

Serious incidents

Whilst serious incidents in relation to medical devices are rare, it is acknowledged that incidents may happen. As an organisation, Cochlear recognises the potential for harm and will respond to any reported serious incident.

What is a serious incident?

A 'serious incident' means any event that directly or indirectly has caused or could have caused an unexpected or unwanted event including any of the following:

- The death of a patient, user or other person,
- The temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
- A serious public health threat.

Reporting a serious incident

There is no definitive list of events/incidents that constitute a serious incident, however all serious incidents should be reported to:

- Your local Cochlear office
www.cochlear.com/intl/contact/global-offices

People within the European Union should also report all serious incidents to:

- Your National Competent Authority
http://ec.europa.eu/growth/sectors/medical-devices/contacts_en

People within Australia should also report all serious incidents to:

- Therapeutic Goods Administration
<https://www.tga.gov.au>

Electromagnetic compatibility (EMC)

Warning

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 in.) to any part of your Nucleus 8 Sound Processor including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Interference may occur in the vicinity of equipment marked with the following symbol:

Environmental protection

Your Cochlear Nucleus 8 Aqua+ coil contains electronic components subject to the Directive 2002/96/EC on waste electrical and electronic equipment.

Help protect the environment by not disposing of your coil with your unsorted household waste. Please recycle or dispose of your coil according to your local regulations for electronics.

Equipment classification

Your sound processor is internally powered equipment Type B applied part as described in the international standard IEC 60601-1:2005/A1:2012, Medical Electrical Equipment–Part 1: General Requirements for Basic Safety and Essential Performance.

Legal statement

The statements made in this version of the User Guide are believed to be true and correct as of the date of publication. However, specifications are subject to change without notice.

© Copyright 2022 Cochlear Limited. All rights reserved.

Trademark legal statement

ACE, Advance Off-Stylet, AOS, Ardiium, AutoNRT, Autosensitivity, Baha, Baha SoftWear, BCDrive, Beam, Bring Back the Beat, Button, Carina, Cochlear, 科利耳, コクレア, 코클리어, Cochlear SoftWear, Contour, コントゥア, Contour Advance, Custom Sound, DermaLock, Freedom, Hear now. And always, Hugfit, Human Design, Hybrid, Invisible Hearing, Kanso, LowPro, MET, MP3000, myCochlear, mySmartSound, NRT, Nucleus, Osia, Outcome Focused Fitting, Off-Stylet, Piezo Power, Profile, Slimline, SmartSound, Softip, SoundArc, True Wireless, the elliptical logo, Vistafix, Whisper, WindShield and Xidium are either trademarks or registered trademarks of the Cochlear group of companies.

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- Administrasi Barang-Barang Terapi (Therapeutic Goods Administration) <https://www.tga.gov.au>

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© Hak Cipta 2022 Cochlear Limited. Hak cipta dilindungi undang-undang.

Pernyataan hukum merek dagang

ACE, Advance Off-Stylet, AOS, Ardium, AutoNRT, Autosensitivity, Baha, Baha SoftWear, BCDrive, Beam, Bring Back the Beat, Button, Carina, Cochlear, 科利耳, コクレア, 코클리어, Cochlear SoftWear, Contour, コントウア, Contour Advance, Custom Sound, DermaLock, Freedom, Hear now. And always, Hugfit, Human Design, Hybrid, Invisible Hearing, Kanso, LowPro, MET, MP3000, myCochlear, mySmartSound, NRT, Nucleus, Osia, Outcome Focused Fitting, Off-Stylet, Piezo Power, Profile, Slimline, SmartSound, Softip, SoundArc, True Wireless, logo elipsnya, Vistafix, Whisper, WindShield dan Xidium adalah merek dagang atau merek dagang terdaftar milik grup perusahaan Cochlear.

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ACE, Advance Off-Stylet, AOS, Ardium, AutoNRT, Autosensitivity, Baha, Baha SoftWear, BCDrive, Beam, Bring Back the Beat, Button, Carina, Cochlear, 科利耳, コクレア, 코클리어, Cochlear SoftWear, Contour, コントウア, Contour Advance, Custom Sound, DermaLock, Freedom, Hear now. And always, Hugfit, Human Design, Hybrid, Invisible Hearing, Kanso, LowPro, MET, MP3000, myCochlear, mySmartSound, NRT, Nucleus, Osia, Outcome Focused Fitting, Off-Stylet, Piezo Power, Profile, Slimline, SmartSound, Softip, SoundArc, True Wireless, 타원형 로고, Vistafix, Whisper, WindShield 및 Xidium은 코클리어 그룹의 상표 또는 등록 상표입니다.

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© Hak cipta 2022 Cochlear Limited. Hak cipta terpelihara.

Pernyataan undang-undang tanda dagangan

ACE, Advance Off-Stylet, AOS, Ardium, AutoNRT, Autosensitivity, Baha, Baha SoftWear, BCDrive, Beam, Bring Back the Beat, Button, Carina, Cochlear, 科利耳, コクレア, 코클리어, Cochlear SoftWear, Contour, コントウア, Contour Advance, Custom Sound, DermaLock, Freedom, Hear now. And always, Hugfit, Human Design, Hybrid, Invisible Hearing, Kanso, LowPro, MET, MP3000, myCochlear, mySmartSound, NRT, Nucleus, Osia, Outcome Focused Fitting, Off-Stylet, Piezo Power, Profile, Slimline, SmartSound, Softip, SoundArc, True Wireless, logo elips, Vistafix, Whisper, WindShield dan Xidium ialah tanda dagangan atau tanda dagangan berdaftar syarikat kumpulan Cochlear.

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ACE, Advance Off-Stylet, AOS, Ardium, AutoNRT, Autosensitivity, Baha, Baha SoftWear, BCDrive, Beam, Bring Back the Beat, Button, Carina, Cochlear, 科利耳, コクレア, 코클리어, Cochlear SoftWear, Contour, コントウア, Contour Advance, Custom Sound, DermaLock, Freedom, Hear now. And always, Hugfit, Human Design, Hybrid, Invisible Hearing, Kanso, LowPro, MET, MicroDrive, MP3000, myCochlear, mySmartSound, NRT, Nucleus, Osia, Outcome Focused Fitting, Off-Stylet, Piezo Power, Profile, Slimline, SmartSound, Softip, SoundArc, True Wireless, the elliptical logo, Vistafix, Whisper, WindShield และ Xidium เป็นเครื่องหมายการค้าหรือไม่มีเครื่องหมายการค้าจดทะเบียนของกลุ่มบริษัท Cochlear

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ACE, Advance Off-Stylet, AOS, Ardium, AutoNRT, Autosensitivity, Baha, Baha SoftWear, BCDrive, Beam, Bring Back the Beat, Button, Carina, Cochlear, 科利耳, コクレア, 코클리어, Cochlear SoftWear, Contour, コントウア, Contour Advance, Custom Sound, DermaLock, Freedom, Hear now. And always, Hugfit, Human Design, Hybrid, Invisible Hearing, Kanso, LowPro, MET, MP3000, myCochlear, mySmartSound, NRT, Nucleus, Osia, Outcome Focused Fitting, Off-Stylet, Piezo Power, Profile, Slimline, SmartSound, Softip, SoundArc, True Wireless, logo hình elip, Vistafix, Whisper, WindShield và Xidium là các nhãn hiệu thương mại hoặc nhãn hiệu thương mại đã đăng ký của tập đoàn Cochlear.

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Cochlear Canada Inc 2500-120 Adelaide
Street West, Toronto, ON M5H 1T1, Canada
Tel: +1 416 972 5082 Fax: +1 416 972 5083

Cochlear Europe Ltd 6 Dashwood Lang
Road, Bourne Business Park, Addlestone,
Surrey KT15 2HJ, United Kingdom
Tel: +44 1932 26 3400 Fax: +44 1932 26 3426

Cochlear Benelux NV Schaliënhoeverdreef
20 i, B-2800 Mechelen, Belgium
Tel: +32 15 79 55 11 Fax: +32 15 79 55 70

Cochlear France S.A.S. 135 Route de
Saint-Simon, 31035 Toulouse, France
Tel: +33 5 34 63 85 85 (International) or 0805
200 016 (National) Fax: +33 5 34 63 85 80

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Cochlear (HK) Limited Room 1404-1406,
14/F, Leighton Centre, 77 Leighton Road,
Causeway Bay, Hong Kong
Tel: +852 2530 5773 Fax: +852 2530 5183

Cochlear Korea Ltd 2nd Floor, Yongsan
Centreville Asterium, 25, Hangang-daero
30 gil, Yongsan-gu, Seoul, Korea (04386)
Tel: +82 2 533 4450 Fax: +82 2 533 8408

Cochlear Medical Device (Beijing) Co., Ltd
Unit 2608-2617, 26th Floor, No.9 Building,
No.91 Jianguo Road, Chaoyang District,
Beijing 100022, P.R. China
Tel: +86 10 5909 7800 Fax: +86 10 5909 7900

**Cochlear Medical Device Company India
Pvt. Ltd** Ground Floor, Platina Building, Plot
No C-59, G-Block, Bandra Kurla Complex,
Bandra (E), Mumbai – 400 051, India
Tel: +91 22 6112 1111 Fax: +91 22 6112 1100

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Cochlear Middle East FZ-LLC Dubai
Healthcare City, Al Razi Building 64, Block A,
Ground Floor, Offices IR1 and IR2, Dubai,
United Arab Emirates
Tel: +971 4 818 4400 Fax: +971 4 361 8925

Cochlear Latinoamérica S.A. International
Business Park, Building 3835, Office 403,
Panama Pacífico, Panama
Tel: +507 830 6220 Fax: +507 830 6218

Cochlear NZ Limited Level 4, Takapuna
Towers, 19-21 Como St, Takapuna, Auckland
0622, New Zealand
Tel: + 64 9 914 1983 Fax: 0800 886 036

P2007044-D2006674 V2

Multilingual version of P1940280-D1891413 V2 2022-04

Cochlear®

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© Copyright 2022 Cochlear Limited. สงวนลิขสิทธิ์