

Cochlear™

Nucleus® CR110 Remote Assistant User Guide



Hear now. And always



Symbols used in this guide



NOTE

Important information or advice.



WARNING (harmful)

Potential safety hazards and serious adverse reactions.

Could cause harm to person.

Contents

Introduction.....	7
Purpose of this guide	7
Cochlear Nucleus CR110 Remote Assistant	7
Getting to know your remote assistant	9
Remote assistant components	9
Getting started.....	11
Charging and recharging your remote assistant battery	11
Charging your remote assistant battery for the first time	11
Recharging your remote assistant battery.....	11
When to recharge your remote assistant battery.....	12
How to recharge your remote assistant battery	13
Turning your remote assistant on and off	18
Locking and unlocking your remote assistant keypad.....	19
Using your remote assistant.....	21
Selecting your remote assistant language.....	21
Pairing your remote assistant with your processor.....	22
Holding and operating your remote assistant	24
Using Simple and Advanced modes.....	26
Simple mode.....	26
Advanced mode	26
Selecting Simple or Advanced mode.....	27
Home screen	28
Selecting programs	29
Increasing and decreasing volume	31
Increasing and decreasing sensitivity.....	32
Fixed volume and sensitivity setting.....	33
Checking the status of your processor	33
Using the telecoil	35
Using Auto Telecoil.....	36
Switching from Auto Telecoil to telecoil.....	36

Turning an audio accessory on and off	36
Understanding remote assistant indicator lights	37
Resetting your processor.....	38
Using the Hard Reset button.....	38
Bilateral control	39
Dual control mode	41
Right or left control modes	41
Switching between control modes.....	42
Adjusting settings	45
Adjusting your processor settings	46
Locking and unlocking your processor buttons	46
Changing the telecoil and microphone mixing ratio	46
Changing the audio accessory and microphone mixing ratio ...	47
Setting your processor indicator light preferences	48
Enabling and disabling private tones	50
Adjusting your remote assistant settings.....	51
Changing your remote assistant speaker volume	51
Enabling and disabling audio and visual alerts	51
Demo mode.....	53
Caring for your remote assistant.....	55
Protecting your remote assistant	55
Using the remote assistant lanyard	56
Cleaning your remote assistant	58
General warnings and precautions for using your remote assistant.....	59
Using your remote assistant to troubleshoot your processor	61
Checking your processor coil.....	62
Other information.....	63
Viewing your remote assistant version	63
Your warranty and registration form.....	63
Technical information.....	63

Specifications	64
Environmental conditions	65
Environmental Protection	66
Equipment classification	66
Certification and applied standards.....	66
FCC (Federal Communications Commission) and Canadian IC compliance.....	67
Labelling symbols.....	68
Legal statement	70

Introduction

Purpose of this guide

This user guide is intended for recipients, parents and carers using the Cochlear™ Nucleus® CR110 Remote Assistant. It provides instructions for using the remote assistant to communicate with the Cochlear™ Nucleus® CP810 Sound Processor.

Cochlear Nucleus CR110 Remote Assistant

Your Cochlear Nucleus CR110 Remote Assistant is a hand-held remote control that enables you to communicate with your Cochlear Nucleus CP810 Sound Processor. It can also receive communication from your processor. If you have two processors, your remote assistant can communicate with both processors. For more information on the CP810 Sound Processor, refer to the CP810 Sound Processor User Guide.

You can use your remote assistant to control the functions and settings of your processor. It provides easy access to key information on the status of your processor (e.g. processor battery status). You can also use your remote assistant to troubleshoot your processor.



Figure 1: CP810 Sound Processor and CR110 Remote Assistant



Figure 2: Remote assistant components

Getting to know your remote assistant

Remote assistant components

1. Indicator light that provides a visual indication of remote assistant function (e.g. remote assistant battery is charging) or a problem (e.g. remote assistant is not communicating with your processor).
2. Upper left soft key to select a listening program or access the processor settings menu.
3. Toggle button to:
 - Turn telecoil (quick press) and Auto Telecoil (press and hold) on and off.
 - Turn audio accessory on and off (quick press).
4. Lower left soft key to select a listening program, or to access the troubleshooting guide.
5. Cochlear key to:
 - Turn the remote assistant on and off (press and hold).
 - Return to the Home screen at any time (quick press).
 - View processor status and menus (quick press).
6. Left directional key to decrease sensitivity (if available) or scroll left.
7. Lower directional key to decrease volume or scroll down.
8. USB cable socket for connecting the remote assistant USB cable.
9. Right directional key to increase sensitivity (if available) or scroll right.
10. Upper directional key to increase volume or scroll up.
11. Lower right soft key to select a listening program, or to access the troubleshooting guide.
12. Remote assistant display (Figure 2 shows Home screen).
13. Remote assistant keypad lock or unlock slider.
14. Upper right soft key to select a listening program or access the remote assistant settings menu.




Figure 3: The back of the remote assistant

- 15. Hole for attaching the remote assistant lanyard.
- 16. Personal identification label.
- 17. Remote assistant speaker.
- 18. Product information label.
- 19. In-built coil sensor guides.
- 20. Hard Reset button to reset the remote assistant if the Cochlear key does not work, or if the remote assistant freezes.

Getting started

This section helps you get started with your new remote assistant. It explains:

- How to charge your remote assistant battery.
- How to turn your remote assistant on and off.
- How to lock and unlock your remote assistant keypad.

Charging and recharging your remote assistant battery

Your remote assistant has a non-removable rechargeable battery.

Charging your remote assistant battery for the first time

Cochlear recommends that you fully charge the battery before using your remote assistant for the first time. It may take approximately 2 to 4 hours to fully charge a completely empty battery. As the battery ages, it may take longer to fully charge the battery.

Recharging your remote assistant battery

A fully charged battery should provide you with approximately 1 to 2 weeks of battery life. Battery life, however, may vary depending on:

- How often and how long you use your remote assistant.
- Whether you turn off your remote assistant when not in use.
- Whether your remote assistant is controlling one or two processors.

When to recharge your remote assistant battery

Charge the battery:

- When the remote assistant battery status screen shows that the battery is low. When you turn on your remote assistant, the battery status screen displays after it displays the start-up (Cochlear logo) screen. The less green bars you see on the battery icon, the less charge you have. The remote assistant keypad lock/unlock screen also indicates battery status. For more information, see *Locking and unlocking your remote assistant keypad*.
- Once a week, even if your remote assistant does not indicate that the battery is low or empty.

It is recommended that you follow these guidelines to maximise the life expectancy of the remote assistant battery.

Recommended temperature range for recharging the battery

To fully charge the remote assistant battery, recharge the battery at a temperature range of 15 °C (59 °F) to 29 °C (84 °F). If the temperature range is lower [between -2.5 °C (28 °F) and 15 °C (59 °F)] or higher [between 29 °C (84 °F) and 52 °C (126 °F)], only 75 to 80 percent of the battery will be charged. If the temperature is too high [over 52 °C (126 °F)] or too low [below -2.5 °C (28 °F)], the battery will not charge.

How to recharge your remote assistant battery

You can recharge the battery using any of the following battery charging equipment available from Cochlear:

- Cochlear™ Nucleus® CR100 Series Global Power Adaptor and Cochlear™ Nucleus® CR100 Series USB Cable.
- Cochlear Nucleus CR100 Series USB Cable attached directly to the USB port on your personal computer.
- Cochlear™ Nucleus® CP800 Series Battery Charger and Cochlear Nucleus CR100 Series USB Cable. The CP800 Series Battery Charger allows you to charge your remote assistant and your processor rechargeable battery module/s at the same time. The Cochlear™ Nucleus® CP800 Series Global Power Adaptor connects the CP800 Series Battery Charger to mains power.



WARNING

Ensure that you only use battery charging equipment supplied by Cochlear to charge the remote assistant battery.

To recharge your remote assistant battery using the CR100 Series Global Power Adaptor and the CR100 Series USB Cable:

1. Connect the CR100 Series USB Cable to the USB socket at the base of the remote assistant.



2. Plug the other end of the CR100 Series USB Cable into the CR100 Series Global Power Adaptor.



3. Connect the CR100 Series Global Power Adaptor to a mains power outlet.

To recharge your remote assistant battery using your computer and the CR100 Series USB Cable:

1. Connect the CR100 Series USB Cable to the USB socket at the base of the remote assistant.



2. Plug the other end of the CR100 Series USB Cable into the USB port of your computer.



When you connect the CR100 Series USB cable to the computer, a 'New hardware found' wizard (dialog) may display on your computer. If you have Internet access, click 'Next' on your computer to install the software for the Cochlear CR110 device.

If you do not have Internet access, press Cancel to dismiss the wizard.

**NOTE**

If your computer is hibernating, in standby mode, or shuts down while charging, the remote assistant battery will stop charging.

To recharge your remote assistant battery using the CP800 Series Battery Charger and the CR100 Series USB Cable:

1. Connect the CR100 Series USB Cable to the USB socket at the base of the remote assistant.



2. Connect the other end of the CR100 Series USB Cable to the CP800 Series Battery Charger.



3. Connect the CP800 Series Battery Charger to the CP800 Series Global Power Adaptor.
4. Connect the CP800 Series Global Power Adaptor to a mains power outlet.

When the CR100 Series USB Cable is first plugged into the CP800 Series Battery Charger, a battery charging animation will appear on your remote assistant display.

When the remote assistant battery is charging, the indicator light on your remote assistant is a steady orange. Once the battery is fully charged, the indicator light is a steady green and the battery icon on the battery status screen displays four green bars.



Figure 4: Battery status screen showing battery is fully charged

Turning your remote assistant on and off

To turn on your remote assistant, press and hold the Cochlear key until the start-up screen (Cochlear logo) appears on the display.



Figure 5: Pressing the Cochlear key



Figure 6: Start-up screen

To turn off your remote assistant, press and hold the Cochlear key until the display turns off after displaying the Goodbye screen.



Figure 7: Goodbye screen

To save battery life, your remote assistant display goes blank after approximately 25 seconds. To 'wake up' the display, you can press any key on your remote assistant.

Locking and unlocking your remote assistant keypad

You can lock your remote assistant keypad to prevent children from changing the settings, or to avoid accidental key presses changing your settings (e.g. when the remote assistant is in your handbag or pocket).

You can lock the keypad using the lock/unlock slider on the remote assistant.

To lock the remote assistant keypad, slide the lock/unlock slider downwards. The lock/unlock status screen displays when the keypad is locked. This screen also shows the remote assistant battery status. You cannot turn off your remote assistant when the keypad is locked.

To unlock, slide the lock/unlock slider upwards.



Figure 8: Locking/unlocking the remote assistant keypad



Figure 9: Remote assistant keypad locked



Figure 10: Remote assistant keypad unlocked

Using your remote assistant

This section describes how to use your remote assistant. It describes:

- How to select your remote assistant language.
- How to pair your remote assistant with your processor.
- How to hold and operate your remote assistant.
- How to access and perform key functions.
- The indicator lights on your remote assistant and what they mean.

Selecting your remote assistant language

When you turn on your remote assistant, the Language Select screen displays if your remote assistant is not paired with your processor. This screen allows you to select the language you wish to use.



Figure 11: Language Select screen

To change the language, scroll left or right using the left and right directional keys. You can also use the lower left and lower right soft keys.

To pair your remote assistant after changing the language, place the coil on the back of your remote assistant. For more information, see the following section.

Pairing your remote assistant with your processor

Pairing allows your remote assistant to communicate with your processor. If you are a bilateral user, you can pair the remote assistant with up to two processors at the same time.

You need to pair your remote assistant if:

- Your remote assistant is not already paired with your processor.
- You have received a new processor.
- You have received a second processor (bilateral user).
- You have received a new remote assistant.
- You have unpaired your remote assistant from your processor.
- Your remote assistant or processor firmware version has been updated.

To pair your remote assistant with your processor:

1. Ensure that your processor and remote assistant are turned on.
2. Place your coil on the back of your remote assistant, ensuring the coil is placed within the coil sensor guides.



Figure 12: Placing the coil on the coil sensor guides

The coil check and pair screen displays.



Figure 13: Coil check and pair screen

3. Press the lower right soft key to pair your remote assistant. On successful pairing, you will see a confirmation screen.

A green tick (✓) on the confirmation screen indicates that the remote assistant has been successfully paired with your processor. If your processor is on the left side, a blue processor icon displays. If your processor is on the right side, a red processor icon displays.

If using two processors (bilateral user), repeat steps 2 and 3 to pair your second processor.

Once the remote assistant is paired, the Home screen displays.



NOTE

If you pair a second processor to the same side, your remote assistant unpairs your first processor and pairs the second processor.

Holding and operating your remote assistant

To optimise communication between your remote assistant and your processor:

- Hold your remote assistant in the palm of your hand in front of you.



Figure 14: Holding the remote assistant in the palm of the hand

- Keep your remote assistant and your processor within two meters (approximately six feet) of each other. Operating it at a greater distance may result in either increased delay between processor and remote assistant response time, or loss of communication.
- Hold or position your remote assistant on the same side of the body as your processor (if communicating with one processor), or between the two processors (if communicating with two processors).
- It is not necessary to point your remote assistant at your processor.

When your remote assistant is within operating distance of your processor and they are communicating, the indicator light on your remote assistant is a steady blue.

You will see blue indicator light flashes and the following icon on the remote assistant display if:

- Your remote assistant is not within operating distance of your processor, or
- There is interference, or
- Your processor is not available (e.g. processor is turned off).



Figure 15: Indicates remote assistant is not within operating distance, there is interference, or processor is not available

Sources of interference could include:

- 2.4 GHz cordless digital telephones.
- Bluetooth devices such as hands free mobile phone kits, wireless keyboards and wireless mice.
- Wireless computing networks such as WiFi.
- Game consoles with 2.4 GHz wireless links.

If you experience interference in the communication between your remote assistant and your processor, try the following:

- Hold your remote assistant in a different position.
- Move away from the vicinity of any electronic equipment that may be causing interference.

If you continue to experience interference, consult your clinician.

Using Simple and Advanced modes

You can operate your remote assistant in Simple or Advanced mode.

Simple mode

Simple mode allows you to do the following:

- Select a language for your remote assistant.
- Pair your remote assistant with your processor.
- Change your listening program.
- Change volume.
- Change microphone sensitivity (if available).
- Turn the telecoil/auto telecoil on and off.
- Change the Home screen display from program names to program numbers.
- View your remote assistant and processor status.
- View visual alerts.
- Reset your processor to the values set by your clinician.
- Use Demo mode.

Advanced mode

Advanced mode allows access to the following additional features on your remote assistant:

- Adjust processor settings (e.g. enable private tones)
- Adjust remote assistant settings (enable visual alerts)
- Processor troubleshooting

Selecting Simple or Advanced mode

Your remote assistant is automatically in Simple mode when you turn on your remote assistant for the first time. You can access Advanced mode, by pressing the upper left soft key when the remote assistant battery status screen displays.



Figure 16: Accessing Advanced mode (1) or Simple mode (2)

To switch between Advanced and Simple modes:

1. Turn off your remote assistant and turn it back on.
2. When the battery status screen displays, press the upper left soft key to select Advanced mode, or the upper right soft key to select Simple mode.

Home screen

The Home screen allows you to select listening programs and perform key functions (e.g. change volume).

Up to four programs are available for your use, depending on the number of programs enabled on your processor. Programs help you hear better in different listening situations.

The available programs can include:

- Everyday program: This is ideal for the typical listening situations of everyday life. 
- Noise program: This is suitable for environments with ongoing background noise coming from all sides. 
- Focus program: This is useful where there is significant background noise, but the focus is on hearing what one person or a small group of people are saying. 
- Music program: This can be used to listen to all types of music from any audio source. 

When you turn on your remote assistant, it takes you to the Home screen after displaying the start-up screen (Cochlear logo) and the remote assistant battery status. The Home screen displays only when your remote assistant is paired with your processor.

Selecting programs

To select a program, press the soft key corresponding to the program you want to use. For example, if you want to use the Noise program, press the lower left soft key (see Figure 17).



Figure 17: Soft keys for selecting a program from Home screen

The Home screen indicates your current program by displaying the icon for that program. For example, if your current program is Everyday, the icon for that program displays.




Figure 18: Home screen displaying program icon for Everyday

To return to the Home screen at any time, press the Cochlear key.

When you turn on the telecoil or connect an audio accessory to your processor, an icon displays next to your current program icon.




Figure 19: Home screen displaying icons for telecoil (1), Auto Telecoil (2) and audio accessory (3)

Switching between program names and program numbers

If you are more familiar with program numbers (e.g. if you have previously used a Freedom™ processor), you can press and hold the upper left and upper right soft keys at the same time from the Home screen to display the program numbers for your programs (e.g. If you only have an Everyday and Noise program enabled on your processor, these two program names will be replaced with program numbers 1 and 2).

To return the display to program names, press and hold the upper left and upper right soft keys at the same time until the names display.

Increasing and decreasing volume

Volume controls your perception of loudness. You can:

- Reduce the volume if the sound from the microphone, telecoil or audio accessory is uncomfortably loud.
- Increase the volume if the sound from the microphone, telecoil or audio accessory, including your own voice, is too soft.

To increase the volume for your current program, press the upper directional key from the Home screen.

To decrease the volume, press the lower directional key from the Home screen.



Figure 20: Directional keys for increasing volume (1) and decreasing volume (2)



Figure 21: Orange tab on volume scale indicates default volume (3)

You can set volume between 1 and 10. The orange tab (3) on the volume scale indicates the default volume set by your clinician. To return to the default volume at any time, change the volume to the level indicated by this tab.

Increasing and decreasing sensitivity



NOTE

Sensitivity may not be available on your sound processor. Please contact your clinician for more information.

Sensitivity controls the softest level of sound picked up by the microphone, telecoil or an audio accessory. You can:

- Reduce the sensitivity of the sound to reduce background noise in noisy situations.
- Increase the sensitivity of the sound to hear very soft sounds in quiet situations.

To increase the sensitivity for your current program, press the right directional key from the Home screen.

To decrease sensitivity, press the left directional key from the Home screen.




Figure 22: Directional keys for decreasing sensitivity (1) and increasing sensitivity (2)




Figure 23: Orange tab on Sensitivity scale indicates default sensitivity

You can set sensitivity between 0 and 20. The orange tab (3) on the sensitivity scale indicates the default sensitivity level set by your clinician. To return to the default sensitivity level at any time, change the sensitivity to the level indicated by this tab.

Fixed volume and sensitivity setting

Your clinician can also set volume and sensitivity at a fixed level so they cannot be changed using your remote assistant.

Checking the status of your processor

Your remote assistant allows you to check if your processor is functioning properly. You can also check the status of the processor battery, and the level of audio input from the microphone, telecoil or an audio accessory. If you are using more than one source of audio input (e.g. mix of audio accessory and microphone input), the remote assistant indicates the level of input from the combined sources.

To view the status of your processor, press the Cochlear key from the Home screen. The processor status screen displays only when your remote assistant is paired with your processor.

If your processor is functioning properly, the status screen displays a green processor icon.



Figure 24: Processor status screen indicating processor is functioning properly

If your remote assistant is receiving any alerts from your processor (e.g. a coil error), a visual alert displays.




Figure 25: Coil off visual alert

If a visual alert displays, you can press one of the lower soft keys to read more information about the problem and how to fix it. The first processor status screen displays for a few seconds before switching to a screen displaying battery status and audio source.

This screen shows the current processor battery status and audio input level. The more green bars you see on the battery icon, the more charge is left in the processor battery. The audio meter shows the current level of audio input. The icon next to the audio meter shows the source of audio input (e.g. microphone).




Figure 26: Processor battery status (1) source of audio input (2) and level of audio input (3)

If your processor is using a Standard Battery Module (zinc air batteries), an hourglass icon displays on the battery status screen for the first ten minutes. This indicates that the processor is checking battery status.



Figure 27: Processor is checking battery status

Using the telecoil

Your processor has an in-built telecoil. The telecoil helps reduce background noise when using the telephone, or when you are in a roomloop or neckloop environment. Your clinician needs to enable the telecoil feature on your processor.

If the telecoil feature is enabled, you can turn on the telecoil by pressing the Toggle button on the remote assistant. The telecoil icon will display next to your current program icon on the Home screen.

To turn off the telecoil, press the Toggle button again. The telecoil icon disappears from the Home screen.



Figure 28: Home screen with telecoil turned on in the Everyday program



Figure 29: Home screen with telecoil turned off in the Everyday program

Using Auto Telecoil

Auto Telecoil is a feature that can be offered by your clinician. When you enable Auto Telecoil, the telecoil automatically turns on when you use the telephone, or when you are in a roomloop environment. The telecoil automatically turns off when:

- You move the telephone away from your processor (e.g. when the telephone call ends)
- You exit the roomloop environment.

It may take up to ten seconds for the telecoil icon to disappear from the screen.

To enable Auto Telecoil, press and hold the Toggle button on the remote assistant until you see the Auto Telecoil icon on the Home screen. When your telephone call starts, the colour of the Auto Telecoil icon changes to black, indicating that you are using the telecoil.



Figure 30: Auto Telecoil turned on (1) and in use (2)

Switching from Auto Telecoil to telecoil

To switch from Auto Telecoil to telecoil, press the Toggle button. The telecoil icon will display on the Home screen. To turn off telecoil, press the Toggle button again.

Turning an audio accessory on and off

To turn on and off an audio accessory connected to your processor, press the Toggle button on the remote assistant.

Understanding remote assistant indicator lights

You may see the following indicator lights on your remote assistant:

Indicator lights	What they mean
 Steady blue	The remote assistant is within operating distance of the processor and the remote assistant and processor are communicating.
 Flashing blue	The processor is not communicating with the remote assistant. For example, they are not within operating distance of each other, or they are not paired.
 Flashing orange	<ul style="list-style-type: none"> • A visual alert needs your attention (e.g. a coil error). • The remote assistant battery is low or empty. • The remote assistant firmware may need an update.
 Steady orange	The remote assistant battery is charging.
 Steady green	The remote assistant battery is fully charged.

Resetting your processor

If you accidentally change a processor setting, or you wish to return to the values set by your clinician, you can reset your processor.

To reset your processor to the values set by your clinician:

1. Press and hold both the lower left and lower right soft keys at the same time from the Home screen. The processor reset screen displays.
2. Press the lower right soft key to confirm that you want to reset your processor.



Figure 31: Processor reset screen

Using the Hard Reset button

If the Cochlear key does not work or your remote assistant display freezes, you can press the Hard Reset button on the back of the remote assistant to reset your remote assistant.




Figure 32: Hard reset button on the remote assistant

Bilateral control

If you have two implants (bilateral user), you can pair your remote assistant with two processors. This allows you to control both processors at the same time.

When your remote assistant is paired with two processors, the remote assistant display is divided vertically in two on the following screens:

- Home screen.
- Status screens.
- Visual alerts screen.




Figure 33: Bilateral mode Home screen



Figure 34: Bilateral mode processor battery status screen




Figure 35: Bilateral mode visual alerts screen

From the Home screen, you can choose one of three bilateral control modes:

- Dual control mode controls both processors.




Figure 36: Bilateral control mode dual left side (1) and right side (2)

- Left control mode* controls only your left processor.




Figure 37: Bilateral control mode left side

- Right control mode* controls only your right processor.




Figure 38: Bilateral control mode right side

* Advanced mode only

Dual control mode

When you pair two processors, your remote assistant is automatically in dual control mode.

In dual control mode, you can control both processors at the same time to:

- Change programs.
- Change volume or sensitivity (if available).
- Turn on the telecoil or Auto Telecoil.
- Connect an audio accessory.

To change these settings, both processors must be on and your remote assistant must be within operating distance of your processors.

Changes to processor settings (e.g. changes to telecoil and microphone mixing ratio) affect both processors, regardless of the control mode.

Right or left control modes

Sometimes you may want to change a setting on one processor only. For example:

- Turn on the telecoil or Auto Telecoil on the left processor only.
- Change the volume or sensitivity (if available) on the right processor only.

Switching between control modes

To switch between control modes, press and hold the left or right directional key.



Figure 39: Switching from dual control mode to left control mode




Figure 40: Switching from left control mode to dual control mode




Figure 41: Switching from dual control mode to right control mode




Figure 42: Switching from right control mode to dual control mode




Figure 43: Switching from left control mode to right control mode




Figure 44: Switching from right control mode to left control mode

Adjusting settings

You can adjust your processor and remote assistant settings to suit your listening needs. You can change the processor settings only when the processor is turned on and within operating distance. If the remote assistant is paired with two processors, both processors need to be turned on and within operating distance.

You can adjust the settings only when operating your remote assistant in Advanced mode.

To adjust your processor or remote assistant settings:

1. From the Home screen, press the Cochlear key. The processor status screen displays.
2. To access the processor settings menu, press the upper left soft key. To access the remote assistant settings menu, press the upper right soft key.



Figure 45: Accessing Processor settings menu (1) or remote assistant settings menu (2)

3. To select the processor or remote assistant setting you wish to change, press the lower directional key until the desired settings screen displays.

Processor settings that you change apply to all programs, regardless of the program that is selected. When the remote assistant is paired with two processors, all processor settings apply to both processors.

Adjusting your processor settings

You can change the following processor settings:

- Locking and unlocking processor buttons.
- telecoil and microphone mixing ratio.
- audio accessory and microphone mixing ratio.
- processor indicator light preferences.
- Enabling and disabling private tones.

Locking and unlocking your processor buttons

You can lock the two buttons on your processor to prevent children from changing the settings, or to avoid accidental button presses changing your processor settings.



Figure 46: Processor buttons unlocked



Figure 47: Processor buttons locked

To lock the buttons, press the lower right soft key from the lock or unlock screen. To unlock, press the lower left soft key.

Changing the telecoil and microphone mixing ratio

Mixing allows you to combine input from more than one audio input source (e.g. telecoil and microphones).

The telecoil and microphone mixing ratio specifies the ratio of the telecoil input in relation to the microphone input. To focus more on the telecoil input and reduce background sounds, increase the mixing ratio. To hear more background sounds while using the telecoil, decrease the mixing ratio.

You can set the mixing ratio for the telecoil and microphone to any of the following:

- **1:1** (50% telecoil input: 50% microphone input).
- **2:1** (67% telecoil input: 33% microphone input).
- **3:1** (75% telecoil input: 25% microphone input).
- **4:1** (80% telecoil input: 20% microphone input).
- **5:1** (83% telecoil input: 17% microphone input).
- **6:1** (86% telecoil input: 14% microphone input).
- **100%** telecoil input.

When you set telecoil input to 100%, you see the letter 'T' on the mixing ratio screen.

To focus more on the telecoil input, press the lower left soft key from the telecoil and microphone mixing screen. To hear more background sounds, press the lower right soft key.



Figure 48: Letter 'T' indicates Telecoil input is set to 100%

Changing the audio accessory and microphone mixing ratio

The audio accessory and microphone mixing ratio specifies the ratio of audio accessory input in relation to microphone input. To focus more on audio accessory input and reduce background sounds, increase the mixing ratio. To hear more background sounds while using audio accessories, decrease the mixing ratio.

You can set the mixing ratio for an audio accessory and the microphone to any of the following:

- **1:1** (50% audio accessory input: 50% microphone input).
- **2:1** (67% audio accessory input: 33% microphone input).
- **3:1** (75% audio accessory input: 25% microphone input).
- **4:1** (80% audio accessory input: 20% microphone input).
- **5:1** (83% audio accessory input: 17% microphone input).
- **6:1** (86% audio accessory input: 14% microphone input).
- **100%** audio accessory input.

When you set audio accessory input to 100%, you see the letter 'A' on the mixing ratio screen.

To focus more on the audio accessory input, press the lower left soft key from the audio accessory and microphone mixing screen. To hear more background sounds, press the lower right soft key.




Figure 49: Letter 'A' indicates audio accessory input is set to 100%

Setting your processor indicator light preferences

The indicator light on your processor provides a visual confirmation when you change a processor setting, or when there is a problem with the processor.

The processor settings menu allows you to set indicator light preferences for your processor. For visual confirmation that the processor is receiving sound, turn on the sound ok indicator light. For visual confirmation of setting changes and processor problems, turn on the alert and key press indicator light.

To set your indicator light preferences:

1. Select the desired indicator light settings screen (e.g. alerts and key press indicator light settings screen).
2. To turn on the indicator light, press the lower right soft key. To turn it off, press the lower left soft key.



Figure 50: Sound ok indicator light is off




Figure 51: Sound ok indicator light is on



Figure 52: Alert and key press indicator light is off




Figure 53: Alert and key press indicator light is on

Enabling and disabling private tones

A private tone is a sound that only the recipient hears when a setting is changed, or when there is a problem with the processor. To hear a tone when a setting is changed, or when there is a problem, enable processor private tones. For silent use, disable processor private tones.

To enable private tones, press the lower right soft key from the private tones screen. To disable private tones, press the lower left soft key.



Figure 54: Private tones disabled




Figure 55: Private tones enabled

Adjusting your remote assistant settings

You can change the following remote assistant settings to suit your needs:

- Remote assistant speaker volume.
- Remote assistant audio and visual alerts.
- Remote assistant language.

Changing your remote assistant speaker volume

Your remote assistant allows you to change the volume of the sounds associated with visual alerts and key presses.

To increase the speaker volume, press the lower right soft key from the speaker volume screen. To decrease the volume, press the lower left soft key.




Figure 56: Remote assistant speaker volume screen

Enabling and disabling audio and visual alerts

When your remote assistant visual alerts are enabled, an alert pops up on your remote assistant display when there is a problem (e.g. when the coil is off). You will also hear a sound.

If your remote assistant alerts are disabled, you will only see visual alerts when you enter the processor status screen. To save battery power, there are no visual alerts when your remote assistant is in standby mode. The orange indicator light will flash when there is a problem, even if visual alerts are disabled.

To enable alerts, press the lower right soft key from the audio and visual alerts screen. To disable alerts, press the lower left soft key.



Figure 57: Remote assistant audio and visual alerts disabled




Figure 58: Remote assistant audio and visual alerts enabled

Demo mode

Demo mode allows you to familiarise yourself with the features on your remote assistant before you start using it to communicate with your processor.

You can familiarise yourself with bilateral (uses two processors) advanced functions or unilateral (uses one processor) simple functions. Demo mode will not display visual alerts, or connect audio accessories.

To enter Demo mode when your remote assistant is already paired with your processor:

1. Press and hold the top right and bottom left soft keys at the same time from the Home screen. The Demo mode screen displays.
2. To select bilateral advanced functions, press the lower left soft key. To select unilateral simple functions, press the lower right soft key.




Figure 59: Access bilateral advanced functions (1) or unilateral simple functions (2)

To enter Demo mode when your remote assistant is not paired, press and hold the top right and bottom left soft keys at the same time from the pair prompt screen. When your remote assistant is not paired, the pair prompt screen displays if you press the Cochlear key.



Figure 60 Pair prompt screen

When using Demo mode, the word 'Demo' will display on all screens. The indicator light on your remote assistant will also flash green, blue and orange.

To stop using Demo mode, turn off your remote assistant and turn it back on.

You will automatically exit Demo mode after 15 minutes of use, or if you start pairing your remote assistant.

Caring for your remote assistant

Protecting your remote assistant

To protect your remote assistant against impact and scratches, Cochlear recommends that you use a protective case. The following protective cases are available from Cochlear:

- Cochlear™ Nucleus® CR100 Series Belt Case: You can insert the remote assistant into the belt case and attach it to your belt. The belt case can be rotated to adjust its position on the belt. When you want to use your remote assistant, remove it from the belt case.



Figure 61: Remote assistant belt case

- Cochlear™ Nucleus® CR100 Series Leather Case: The leather case provides you with the flexibility of either attaching the remote assistant to your belt, or carrying it in your handbag or pocket. You can use and charge the remote assistant while it is in the leather case.



Figure 62: Remote assistant leather case

- Cochlear™ Nucleus® CR100 Series Protective Case: This is a flexible rubber case into which you can insert the remote assistant. It is suitable for use by both children and adults, and is available in a choice of colours. The protective case can be worn around the neck using the remote assistant lanyard. You can use and charge the remote assistant while it is in the protective case.



Figure 63: Protective case



NOTE

When using the remote assistant while it is in a protective case, you may notice a slight decrease in the operating range. If this happens, either remove the remote assistant from the case when using it, or hold the case close to your processor.

When the remote assistant is attached to your belt, certain arm positions (e.g. arm swinging past the belt when walking) may affect communication. If this happens, remove the remote assistant from the case when using it.

Using the remote assistant lanyard

The Cochlear™ Nucleus® CR100 Series Lanyard allows you to wear the remote assistant around your neck.

You can attach the lanyard to the remote assistant using the hole located on the back of the remote assistant. When threading the lanyard, always insert the lanyard loop through the bottom hole, which is the hole just above the personal identification label.

Thread the lanyard as shown below:

1. Hold the remote assistant in one hand.
2. Insert the lanyard loop through the bottom hole with the thumb and index finger of the other hand.
3. Insert the lanyard strap through the lanyard loop as shown.
4. Pull the lanyard strap upwards until the loop tightens into a knot.



NOTE

Components described in this section may not be available in all countries. Please contact your clinician for more information.

Cleaning your remote assistant

To keep your remote assistant clean, wipe its display and buttons with a soft dry cloth. Do not use alcohol, household cleaners or abrasive materials to clean your remote assistant.



Figure 64: Wiping the remote assistant with a soft dry cloth

If your belt case, leather case or protective case gets dirty, wipe it with a soft dry cloth.

If your remote assistant gets wet:

1. Turn your remote assistant off.
2. Wipe it with a soft dry cloth.
3. Dry your remote assistant for approximately 12 hours. It is not recommended to use a dry aid kit to dry your remote assistant. Do not use an oven or microwave to dry your remote assistant.

To check if your remote assistant is still functioning, turn it on again after it is dry. If your remote assistant no longer functions, return it to your clinician. Cochlear cannot guarantee that it can repair a water-damaged remote assistant.

General warnings and precautions for using your remote assistant

- No modification of this equipment is allowed.
- Parents and carers are advised that unsupervised use of long cables (e.g. coil or accessory cables) may present a risk of strangulation.
- Using your remote assistant cables in any way contradictory to their intended purpose (e.g. wearing the USB cable around the neck) can cause injury.
- When operating machinery, ensure that your remote assistant lanyard does not get caught in the machine.
- If your remote assistant becomes unusually warm, stop using it immediately and seek advice from your clinician.
- Do not operate your remote assistant while driving.
- Do not leave your remote assistant in very hot areas (e.g. in direct sunlight, behind a window, in a car, etc.).
- Avoid placing your remote assistant in environments where excessive moisture and dirt can damage it (e.g. swimming pools).
- Be careful when placing your remote assistant in a bag or pocket with other objects as the display may scratch. You may want to use commercially available protective films to protect the remote assistant display from scratches.
- When charging your remote assistant from mains power, only use the CR100 Series Global Power Adaptor, or the CP800 Series Global Power Adaptor (if using the CP800 Series Battery Charger).
- Do not use excessive force when connecting your remote assistant to the remote assistant USB cable.
- Avoid applying unnecessary pressure to the display, keypad or back of your remote assistant.
- Avoid bending your remote assistant in any way.
- Do not attempt to open or dismantle your remote assistant. If the remote assistant is not working, return it to your clinician.
- When your remote assistant is not in use, turn it off and store it safely.

- Dispose of electrical components in accordance with your local regulations.
- When charging the remote assistant with a computer, ensure the computer conforms to the IEC 60950 standard.

Using your remote assistant to troubleshoot your processor

The following features on your remote assistant can help you troubleshoot your processor:

- Visual alerts that appear on your remote assistant display.




Figure 65: Coil off visual alert

- Your processor status indicator.
- The processor troubleshooting guide on your remote assistant when operating the remote assistant in Advanced mode. You can access the troubleshooting guide by pressing one of the lower soft keys from the processor status screen.

For information on how to troubleshoot your processor, refer to the CP810 and CR110 Troubleshooting Guide provided with your processor documentation.

Checking your processor coil

Your remote assistant has an in-built coil sensor that allows you to check if the coil connected to your processor is receiving sound signals from your processor.

To check the coil using the coil sensor:

1. Remove the coil from the head. Do not use the coil sensor when the coil is on the head.
2. Place the coil on the coil sensor guides on the back of the remote assistant.

If the coil is receiving sound signals from your processor, the coil check and pair screen displays. This only indicates that the coil is receiving sound signals and does not confirm that the coil is functioning properly.

If the coil is not receiving sound signals, the coil check and pair screen does not display.

Once you move the coil away from your remote assistant, the remote assistant returns to its previous display.

Other information

Viewing your remote assistant version

The version of your remote assistant is displayed on the bottom of the start-up screen (Cochlear logo) when you turn your remote assistant on.



Figure 66: Start-up screen showing remote assistant version (1)

Your warranty and registration form

Your warranty is included in the document pack you received with your processor.

Please complete the registration form and return to Cochlear within 30 days of receiving the product.

Technical information

The remote assistant wireless communication link operates in the 2.4 GHz ISM band using GFSK (Gaussian Frequency Shift Keying) on 10 channels. The link uses a proprietary bi-directional communication protocol and operates over a distance of up to 2m from the CP810 processor. When interference is present, the wireless communication link switches between the 10 channels to find a channel where the interference least affects the operation of the link. The remote assistant indicates via its display and light indicator when the processor is not within operating distance, and when the link has been interrupted due to interference. For more information, see *Holding and operating your remote assistant* in this user guide.

Battery life, charge cycles and lifespan

- Battery life means the time a device will run before the disposable batteries must be replaced, or the rechargeable batteries recharged.
- Battery charge cycle means a cycle of one full charge and discharge of the rechargeable battery.
- Battery lifespan means the total number of charge cycles a rechargeable battery will last before the battery life degrades to 80% of its original fully-charged capacity.

Specifications

Physical configuration

The remote assistant comprises:

- Custom analogue and digital integrated circuits based on a micro-processor design, with bi-directional wireless communication capabilities.
- Push-buttons/keys and a slider for control of processor features and its own system features.
- A display for indicating system status and activity.
- A speaker for sounding an audible alert.
- In-built coil sensor to check if the processor's coil is receiving sound signals from the processor.
- An indicator light for visual indication of system function or alert.
- On-board battery charging circuit.
- A temperature sensor for battery safety.
- A mini type B USB connector for charging the device.

Materials

The remote assistant housing, including non-metallic buttons/keys, is made of polycarbonate. The spine and metallic buttons/keys are made of chrome plated stainless steel.

Dimensions

Length	Width	Depth
110 mm	45 mm	12 mm

Weight

The remote assistant weighs approximately 64 g.

Operating characteristics

Characteristic	Range
Wireless technology	Proprietary Low Power Bidirectional Wireless Link
RF frequency	2.4 GHz
Operating battery voltage	3 V to 4.2 V
Battery charging voltage	4.5 V to 5.5 V
Battery capacity	680 mAh (Internal Lithium-ion battery)
Push-key functions	Remote assistant on and off, telecoil and Auto Telecoil on and off, audio accessory on and off, program change, processor and remote assistant settings change
Wireless transmission range	Up to 2 m
Display	128 x 128 pixels, 1.5 inch backlit colour LCD

Environmental conditions

Condition	Minimum	Maximum
Storage & transport temperature	-10 °C (14 °F)	+55 °C (131 °F)
Operating temperature	+5 °C (+41 °F)	+40 °C (+104 °F)
Operating relative humidity	0% RH	90% RH
Operating pressure	700 hPa	1060 hPa

Environmental Protection

Your remote assistant contains electronic components subject to the Directive 2002/96/EC on waste electrical and electronic equipment. Help protect the environment by not disposing of your remote assistant with your unsorted household waste. Please recycle your remote assistant according to your local regulations.

Equipment classification

Your remote control is an accessory to a medical device as described in the international standard IEC 60601-1:2005, Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance.

Certification and applied standards

Your device fulfills the essential requirements listed in Annex 1 of the EC directive 90/385/EEC on Active Implantable Medical Devices as per the conformity assessment procedure in Annex 2.



The year in which authorisation to affix the CE mark was granted was 2015.

FCC (Federal Communications Commission) and Canadian IC compliance

This device complies with part 15 of the FCC Rules and with RSS-210 of Industry Canada. Operation is subject to the following two conditions:

- This device may not cause harmful interference.
- This device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications made to this equipment not expressly approved by Cochlear Limited may void the FCC authorization to operate this equipment.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

FCC ID number: WTOCR11000

IC ID number: 8039A-CR11000

Labelling symbols

The following symbols may appear on your processor or remote components and/or packaging:

Symbol	Meaning
	Refer to instruction manual
	Specific warnings or precautions associated with the device, which are not otherwise found on the label
	Temperature limits
	Humidity Limitation
	CE registration mark with notified body number
	Authorised representative in the European Community
	Keep dry
	Dispose of electrical components in accordance with your local regulations
 R006WWC0198	Radio compliance certification for Japan
	Radio compliance certification for Australia and New Zealand
	Manufacturer
	Catalogue number
	Serial number
	Batch code
	Date of manufacture

Symbol	Meaning
Rx Only	By prescription
	Recyclable material

Legal statement

The statements made in this guide are believed to be true and correct as of the date of publication. However, specifications are subject to change without notice.

© Cochlear Limited 2017

Cochlear™



Cochlear Ltd (ABN 96 002 618 073) 1 University Avenue, Macquarie University, NSW 2109, Australia
Tel: +61 2 9428 6555 Fax: +61 2 9428 6352

Cochlear Ltd (ABN 96 002 618 073) 14 Mars Road, Lane Cove, NSW 2066, Australia
Tel: +61 2 9428 6555 Fax: +61 2 9428 6352

Cochlear Deutschland GmbH & Co. KG Karl-Wiechert-Allee 76A, 30625 Hannover, Germany
Tel: +49 511 542 770 Fax: +49 511 542 7770

Cochlear Americas 13059 E Peakview Avenue, Centennial, CO 80111, USA
Tel: +1 303 790 9010 Fax: +1 303 792 9025

Cochlear Canada Inc 2500-120 Adelaide Street West, Toronto, ON M5H 1T1, Canada
Tel: +1 416 972 5082 Fax: +1 416 972 5083

Cochlear AG EMEA Headquarters, Peter Merian-Weg 4, 4052 Basel, Switzerland
Tel: +41 61 205 8204 Fax: +41 61 205 8205

Cochlear Europe Ltd 6 Dashwood Lang Road, Bourne Business Park, Addlestone, Surrey KT15 2HJ, United Kingdom
Tel: +44 1932 26 3400 Fax: +44 1932 26 3426

Cochlear Benelux NV Schaliënhoevredreef 20 i, B-2800 Mechelen, Belgium
Tel: +32 15 79 55 11 Fax: +32 15 79 55 70

Cochlear France S.A.S. 135 Route de Saint-Simon, 31035 Toulouse, France
Tel: +33 5 34 63 85 85 (International) or 0805 200 016 (National) Fax: +33 5 34 63 85 80

Cochlear Italia S.r.l. Via Larga 33, 40138 Bologna, Italy
Tel: +39 051 601 53 11 Fax: +39 051 39 20 62

Cochlear Nordic AB Konstruktionsvägen 14, 435 33 Mölnlycke, Sweden
Tel: +46 31 335 14 61 Fax: +46 31 335 14 60

Cochlear Tıbbi Cihazlar ve Sağlık Hizmetleri Ltd. Şti.
Çubuklu Mah. Boğaziçi Cad., Boğaziçi Plaza No: 6/1, Kavacık, TR-34805 Beykoz-Istanbul, Turkey
Tel: +90 216 538 5900 Fax: +90 216 538 5919

Cochlear (HK) Limited Room 1204, 12/F, CRE Building, No 303 Hennessy Road, Wanchai, Hong Kong SAR
Tel: +852 2530 5773 Fax: +852 2530 5183

Cochlear Korea Ltd 1st floor, Cheongwon Building 33, Teheran-ro 8 gil, Gangnam-gu, Seoul, Korea
Tel: +82 2 533 4450 Fax: +82 2 533 8408

Cochlear Medical Device (Beijing) Co., Ltd
Unit 2208-2212, Tower B, Gemdale Building, 91 Jianguo Road, Chaoyang District, Beijing 100022, P.R. China
Tel: +86 10 5909 7800 Fax: +86 10 5909 7900

Cochlear Medical Device Company India Pvt. Ltd.
Ground Floor, Platina Building, Plot No C-59, G-Block, Bandra Kurla Complex, Bandra (E), Mumbai – 400 051, India
Tel: +91 22 6112 1111 Fax: +91 22 6112 1100

株式会社日本コクレア (Nihon Cochlear Co Ltd) 〒113-0033 東京都文京区本郷2-3-7 お茶の水元町ビル
Tel: +81 3 3817 0241 Fax: +81 3 3817 0245

Cochlear Middle East FZ-LLC
Dubai Healthcare City, Al Razi Building 64, Block A, Ground Floor, Offices IR1 and IR2, Dubai, United Arab Emirates
Tel: +971 4 818 4400 Fax: +971 4 361 8925

Cochlear Latinoamérica S.A.
International Business Park, Building 3835, Office 403, Panama Pacifico, Panama
Tel: +507 830 6220 Fax: +507 830 6218

Cochlear NZ Limited
Level 4, Takapuna Towers, 19-21 Como St, Takapuna, Auckland 0622, New Zealand
Tel: + 64 9 914 1983 Fax: 0800 886 036

www.cochlear.com

ACE, Advance Off-Stylet, AOS, AutoNRT, Autosensitivity, Beam, Button, Carina, Cochlear, Cochlear SoftWear, コクレア, Codacs, Contour, Contour Advance, Custom Sound, ESPrint, Freedom, Hear now. And always, Hugfit, Hybrid, inHear, Invisible Hearing, Kanso, MET, MicroDrive, MP3000, myCochlear, mySmartSound, NRT, Nucleus, 科利耳, Off-Stylet, SmartSound, Softip, SPrint, True Wireless, the elliptical logo and Whisper are either trademarks or registered trademarks of Cochlear Limited. Ardium, Baha, Baha SoftWear, BCDrive, DermaLock, EveryWear, Vistafix and WindShield are either trademarks or registered trademarks of Cochlear Bone Anchored Solutions AB.

© Cochlear Limited 2017

Hear now. And always

471830 ISS3 JUN17