Cochlear™ Baha® 7 Sound Processor

User manual part A

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
Sound Processor

User manual part A

This guide is intended for recipients and caregivers using the Cochlear™ Baha® 7 Sound Processor as part of the Cochlear Baha System.

Contents

Introduction	5	Care	16
Overview.....	5	Care and maintenance.....	16
Intended purpose.....	5	IP classification	17
Indications	5		
Contraindications	6	Troubleshooting	17
Intended patient population.....	6	Processor will not turn on	18
Intended clinical benefits	6	Sound is too quiet or muffled.....	18
Intended users	6	Sound is too loud or uncomfortable ..	18
Warranty.....	6	You experience feedback (whistling) ..	18
Use	6	Other information	19
Turn on and off	6	Sound processor and parts.....	19
Sound processor indicators	7	Serious incidents.....	20
Change programs	7	Performance characteristics	20
Adjust volume	8	Environmental conditions	21
Share the experience	8	Environmental protection	21
Power	8	Magnetic Resonance Imaging (MRI) ..	21
Battery type	8	Electromagnetic compatibility (EMC) ..	21
Low battery indication	8	Regulatory information	22
Change the battery.....	9	Equipment classification and	
Tamper-resistant battery door	9	compliance.....	22
		Certification and applied standards...	24
Wear	10	Key to symbols	24
Safety line.....	10		
Flight mode	10	Trademark legal notice	26
For users with two sound processors ..	11		
Wireless devices and streaming.....	11		
Streaming via a mobile phone or			
tablet.....	13		
Audio and visual indicators	14		
Wireless audio and visual signals.....	15		
Paediatric mode	15		

Introduction

Congratulations on your choice of the Cochlear™ Baha® 7 Sound Processor. This manual is full of tips and advice on how to best use and care for your Baha sound processor. Be sure to discuss any questions or concerns that you may have regarding your hearing or use of this system with your hearing care professional.

Overview

Note

Additional illustrations, figures 1–9, can be found on the inside of the cover of this user manual.

Intended purpose

The Cochlear Baha System uses bone conduction to transmit sounds to the cochlea (inner ear) with the purpose of enhancing hearing. The Baha 7 Sound Processor is intended to be used as part of the Cochlear Baha System to pick up surrounding sound and transfer it to the skull bone via a Baha Implant, Baha SoundBand, Baha Softband or Baha SoundArc and can be used unilaterally or bilaterally.

Indications

The Baha System is indicated for patients with conductive hearing loss, mixed hearing loss and single-sided sensorineural deafness (SSD). The Baha 7 Sound Processor is indicated for patients with up to 55 dB HL sensorineural hearing loss (SNHL).

Contraindications

There are no existing contraindications for the Baha 7 Sound Processor.

Intended patient population

The Baha 7 Sound Processor is intended for use by recipients of any age, no health conditions specified.

Intended clinical benefits

Most recipients of a bone conduction hearing solution will experience improved hearing performance and quality of life compared to unaided listening.

Intended users

There are two intended user groups of the Baha 7 Sound Processor;

- 1) The recipient or recipient's carer – the carer could be a parent, nurse, caregiver etc.
- 2) The clinician – the clinician is an audiologist with experience of fitting hearing solutions.

Warranty

The warranty does not cover defects or damage arising from, associated with, or related to the use of this product with any non-Cochlear processing unit and/or any non-Cochlear implant. See the *"Cochlear Baha Global Limited Warranty card"* for more details.

Use

Turn on and off

See *figure 1*

The battery door is used to turn the sound processor on and off.

1. To turn on your sound processor, close the battery door completely.
2. To turn off your sound processor, gently open the battery door until you feel the first "click".

When your sound processor is turned off and back on again, it will return to Program 1 and default volume level. If enabled, audio and/or visual signals will let you know that the device is starting up. See *"Audio and visual indicators"* on page 14.

Sound processor indicators

See *figure 2*

Audio signals and the visual indicator will alert you of changes to your sound processor. For a complete overview see *"Audio and visual indicators" on page 14*.

Change programs

See *figure 3*

You can choose between programs to change the way your sound processor deals with sound. You and your hearing care professional will have selected up to four pre-set programs for your sound processor.

Program 1 _____
 Program 2 _____
 Program 3 _____
 Program 4 _____

These programs are suitable for different listening situations. Ask your hearing care professional to fill in your specific programs on the lines.

1. To change the program, press and release the control button located on the top of your sound processor once.
2. If enabled, audio and visual signals will let you know which program you have changed to. See *"Audio and visual indicators" on page 14*.
3. To change to any of the other programs pre-set by your clinician, repeat above steps until you get confirmation that you are in the desired program.

Note

If you are a bilateral recipient, program changes you make to one device will automatically apply to the second device. This function can be enabled or disabled by your hearing care professional.

Adjust volume

Your hearing care professional has set the volume level for your sound processor.

Note

You can change the program and adjust the volume using the optional Remote Control, Cochlear Baha Remote Control 2, Cochlear Wireless Phone Clip, Baha Smart App or from your compatible smartphone or smart device. See *"Wireless devices and streaming"* on page 11.

Share the experience

See *figure 4*

Family members and friends can "share the experience" of bone conduction hearing using the Cochlear test rod, provided with the sound processor.

1. Turn on your sound processor and attach it on the test rod by tilting it into place. You will feel the snap coupling "click" into the notch on the test rod.
2. Hold the test rod against the skull bone behind an ear. (Ensure you are holding the test rod, and not the sound processor). Plug both ears and listen.

Power

Battery type

The Baha 7 Sound Processor uses a 312 size type hearing aid battery (1.45 Volt zinc air, non-rechargeable). Batteries should be replaced as needed, just as you would with many other electronic devices. Battery life will vary with e.g. daily use, volume levels, wireless streaming, sound environment, program setting, and battery strength.

Low battery indication

If activated, the visual and audio signals will alert you when there is approximately one hour of battery power remaining (at this time you may experience lower amplification). If the battery runs down completely, the sound processor will stop working.

Change the battery

See figure 5

1. To replace the battery, remove your sound processor from the head and hold the sound processor with the front facing down.
2. Gently open the battery door until it is completely open.
3. Remove the old battery and dispose of it according to local regulations.
4. Remove the new battery from the packet and peel away the sticker on the + side.
5. Insert the battery into the battery compartment with the + side facing up.
6. Gently close the battery door.

Warning

Batteries can be harmful if swallowed, put in the nose or in the ear. Be sure to keep your batteries out of reach of small children and other recipients in need of supervision. Before use, verify that the tamper-resistant battery door is properly closed. In the event a battery is accidentally swallowed, or stuck in the nose or in the ear, seek immediate medical attention at the nearest emergency centre.

Note

- To maximise battery life, switch off the sound processor when it is not in use.
- Battery life decreases as soon as the battery is exposed to air (when the plastic strip is removed), so be sure to only remove the plastic strip directly prior to use.
- If a battery leaks, replace it immediately.

Tamper-resistant battery door

See figure 6

To prevent the accidental opening of the battery door, an optional tamper-resistant battery door is available. This is particularly useful to prevent children, and other recipients in need of supervision, from accidentally accessing the battery. Contact your hearing care professional for a tamper-resistant battery door.

To use the tamper-resistant battery door:

1. To unlock and turn off the device, carefully insert the tamper resistant tool or the tip of a pen in the small hole on the battery door and gently open the door.
2. To lock and turn on the device, gently close the battery door until it is completely closed.

Wear

Safety line

See *figure 7*

The safety line is designed to reduce the risk of dropping or losing your processor. You can attach a safety line that clips onto your clothing:

1. Pinch the loop on the end of the safety line between your finger and thumb.
2. Pass the loop through the attachment hole in the sound processor from front to back.
3. Pass the clip through the loop and pull the line tight. Attach the clip to your clothing.

Note

Cochlear recommends connecting the safety line when engaging in physical activities. Children should use the safety line at all times.

Flight mode

See *figure 8*

Activate flight mode in situations when you need to deactivate radio signals (wireless functionality), such as when boarding a flight or other areas where radio frequency emission is prohibited.

To activate flight mode:

1. Open and close the battery door on your sound processor three times (open-close, open-close, open-close) within a 10-second period.
2. If enabled, audio and visual signals will confirm that flight mode is activated. See *"Audio and visual indicators"* on page 14.

Status/action	Audio signal	Visual signal
Start up in flight mode	 10x dual beeps	 4 x dual flashes

Follow these steps to deactivate flight mode:

1. Make sure your sound processor has been running for at least 15 seconds before you attempt to turn off flight mode.
2. To turn off flight mode, open and close the battery door once on your sound processor.
3. Let the sound processor run for another 15 seconds or more before turning it off to be certain that flight mode is deactivated.

For users with two sound processors

To make identification easier, ask your hearing care professional to mark your left and right sound processor with the coloured stickers provided (red for right, blue for left).

Wireless devices and streaming

You can use wireless devices to enhance your listening experience. To learn more about the options available, ask your hearing care professional or visit www.cochlear.com.

Compatible devices are:

- Multi-Mic+
- TV-Streamer+
- Remote Control
- Cochlear Baha Remote Control 2
- Cochlear Wireless Mini Microphone 2+
- Cochlear Wireless TV Streamer
- Cochlear Wireless Phone Clip

The sound processor is compatible with the Remote Control, which lets you change program, activate streaming and change volume without the use of a smartphone. For more information, see the Remote Control user guide.

Pair the sound processor with Multi-Mic+ or TV-Streamer+

1. Activate pairing mode on the wireless device (see the user guide for the wireless device).
2. Hold the sound processor close to the wireless device.
3. An indicator light on the device and an audio signal from the sound processor will confirm successful pairing. See *"Audio and visual indicators"* on page 14.

Pair the sound processor with other wireless devices

1. Activate pairing mode on the wireless device (see the user guide for the wireless device).
2. Hold the sound processor close to the wireless device.
3. Turn off the sound processor by opening the battery door.
4. Turn on the sound processor by closing the battery door. An audio signal from the sound processor will confirm successful pairing.

Status/action	Audio signal	Visual signal
Confirmation wireless device paring	 Ripple tone in upward melody	N/A

Activate streaming for wireless devices

When the sound processor is paired with a wireless device, press and hold the control button on your sound processor until you hear an audio signal.

Status/action	Audio signal	Visual signal
Wireless streaming activated or change from one wireless device to another	 Ripple tone upward melody	 1 long flash followed by 1 short flash

End wireless audio streaming

Press and release (short press) the control button on your sound processor. The sound processor will return to the previously used program.

Toggle between wireless devices

If your sound processor is paired with more than one wireless device, you can toggle between the devices in the different channels by pressing the control button (long press) on your sound processor one, two or three times, until you have selected the device you want.

Note

For additional guidance regarding for example pairing, please refer to the user guide of the relevant wireless device.

Streaming via a mobile phone or tablet

The sound processor has built-in technology for streaming audio directly from a compatible iOS or Android device. Make sure that Bluetooth is activated to stream audio directly from your device. For compatibility details and more information, visit www.cochlear.com/compatibility.

Made for iPhone (MFi)

Your sound processor is a Made for iPhone (MFi) hearing device. This allows you to control your sound processor and stream audio directly from your Apple® devices.

Android streaming

Your sound processor is compatible with Bluetooth LE (Low Energy) Audio and the ASHA (Audio Streaming for Hearing Aids) protocol. This allows you to use the direct audio streaming functions of compatible Android devices.

Baha Smart App

Your sound processor is compatible with the Baha Smart App. You can use the Smart App to control and monitor your sound processor and connected wireless devices. The Baha Smart App is also your gateway to Cochlear's Remote Care services, such as Remote Assist. Download the app from App Store or Google Play and follow the on-screen instructions.

Audio and visual indicators

Your hearing care professional can set up your sound processor to show the following audio and visual signals.

Status/action	Audio signal	Visual signal	Comment
Start up	 5 beeps	 4 seconds steady light	Your hearing care professional can set up the audio signal to be 1, 5 or 10 beeps.
Start up in flight mode	 10x dual beeps	 4 x dual flashes	
Change program	 1-4 beeps	 1-4 flashes	The number of flashes and beeps indicates the number of the current program.
Volume up/down	 1 beep	 1 short flash	
Maximum/ minimum volume	 1 long beep	 1 long flash	

Status/action	Audio signal	Visual signal	Comment
Low battery indication	 2 x 4 beeps	 Repeated series of rapid flashes	

Wireless audio and visual signals

Status/action	Audio signal	Visual signal	Comment
Wireless streaming activated or change from one wireless device to another	 Ripple tone upward melody	 1 long flash followed by 1 short flash	
Confirmation wireless device pairing	 Ripple tone in upward melody	N/A	

Paediatric mode

This optional continuous mode is primarily intended for parents and carers who want to receive a visual feedback from their child's sound processor. It can be activated by your hearing care professional. As the child gets older the mode can also be switched off by your hearing care professional.

Status/action	Visual signal	Comment
Low battery indication	•••••••• •••••••• •••••••• Repeated series of rapid flashes	Continuously repeated or repeated with small pauses.
Flight mode	•• •• •• •• 4 x dual flashes	
Program 1–4	• •• ••• •••• 1–4 flashes depending on the chosen program	
Streaming active	▬ • 1 long flash followed by 1 short flash	

Care

Care and maintenance

Your sound processor is a delicate electronic device. Follow these guidelines to keep it in proper working order:

- For cleaning your sound processor and snap coupling, remove the sound processor from your head and use the Baha sound processor cleaning kit and accompanying instructions. The kit is provided by Cochlear in the sound processor box.
- After exercise, wipe your processor with a soft cloth to remove sweat or dirt.
- If the sound processor gets wet or is exposed to a very humid environment, dry it with a soft cloth, remove the battery and let the processor dry out before inserting a new one.
- Remove your sound processor before applying any hair conditioners, mosquito repellent or similar products.
- Turn off and store the sound processor away from dust and dirt.
- A storage case is provided by Cochlear in the sound processor box.
- Avoid exposing your sound processor to extreme temperatures.
- For long-term storage, remove the battery.

 Caution

Do not use other cleaning methods than recommended by Cochlear.

IP classification

The electronics compartment in your sound processor is protected against damage by dust and by immersion in water. Without the battery, the sound processor was tested for immersion in water for 35 minutes at 1.1 meters depth and achieved an IP68 rating. This means that if you, for example, accidentally drop your sound processor in water, the electronics in the device are protected against malfunction due to water ingress. However, your sound processor has a battery that requires air to operate and malfunctions if wet. The sound processor with battery achieves an IP42 rating. This means there is a possibility that if you, for example, are out in rain or in other humid environments, water can block air supply to the battery causing a temporary malfunction. To avoid temporary malfunction, avoid exposing the sound processor to water and always remove it before swimming or bathing.

If your sound processor becomes wet and malfunctions:

1. Remove your sound processor from the head.
2. Open the battery door and remove the battery.
3. Put your sound processor in a container with drying capsules such as a Dri-Aid Kit, etc. Let your sound processor dry out before inserting a new battery. Drying kits are available from most hearing care professionals.

Troubleshooting

Contact your hearing care professional if you have any concerns regarding the operation or safety of your sound processor, or if the solutions below do not resolve your issue.

Processor will not turn on

1. Try turning the sound processor on again. See *"Turn on and off"* on page 6.
2. Replace the battery. See *"Change the battery"* on page 9.
3. The battery requires air to operate. Ensure that the battery air inlet and/or the battery air holes are not covered.
4. Try a different program. See *"Change programs"* on page 7.

Sound is too quiet or muffled

1. Try turning up the volume using a compatible smartphone or wireless device.
2. Check that the sound processor is not wet. If it is wet, let the sound processor dry before use. See *"Care and maintenance"* on page 16.

Sound is too loud or uncomfortable

1. Try turning down the volume of your sound processor. See *"Adjust volume"* on page 8.

You experience feedback (whistling)

1. Check to ensure that the sound processor is not in contact with items such as glasses or a hat, or in contact with your head or ear. See *figure 9*.
2. Try turning down the volume of your sound processor. See *"Adjust volume"* on page 8.
3. Check that there is no external damage to the sound processor.
4. Check that there is no dirt in the connection to your sound processor.

Other information

Sound processor and parts

- The sound processor is suited for use in a home healthcare environment. The home healthcare environment includes locations such as homes, schools, churches, restaurants, hotels, cars, and airplanes, where equipment and systems are less likely to be administered by healthcare professionals.
- A sound processor will not restore normal hearing and will not prevent or improve a hearing impairment resulting from organic conditions.
- Infrequent use of a sound processor may not enable a recipient to attain full benefit from it.
- The use of a sound processor is only part of hearing rehabilitation and may need to be supplemented by auditory and lip reading training.
- The sound processor is a digital, electrical, medical instrument designed for specific use. As such, due care and attention must be exercised by the recipient at all times.
- A discharge of static electricity can damage the electrical components of the sound processor or corrupt the program in the sound processor. If static electricity is present (e.g. when putting on or removing clothes over the head or getting out of a vehicle), you should touch something conductive (e.g. a metal door handle) before your sound processor contacts any object or person. Prior to engaging in activities that create extreme electrostatic discharge, such as playing on plastic slides, the sound processor should be removed.
- If disruptions keep occurring, please contact your clinician to resolve the issue.
- For wireless functionality, only use Cochlear Wireless devices or compatible smart devices.
- No modification of this equipment is allowed.
- Adult supervision is recommended when the recipient is a child.
- Avoid exposing your sound processor to X-ray radiation.

Warning

The sound processor and removable parts of the system (batteries, battery door, safety line) can be lost or may be a choking or strangulation hazard. Keep out of reach of small children and other recipients in need of supervision.

Warning

Do not use damaged product. A damaged product can expose sharp edges or rough surfaces.

Serious incidents

Any serious incident in relation to your device should be reported to your Cochlear representative and to the medical device authority in your country, if available.

Performance characteristics

Measurements according to IEC 60118-9:2019 2nd Ed.

Output vibratory force level (re. 1 μ N)	Max 121 dB
--	------------

(90 dB SPL input, FOG)	HFA 113 dB
------------------------	------------

Acousto-mechanical sensitivity level (re. 1 μ N /	Max 48 dB
---	-----------

20 μ Pa) (50 dB SPL input, FOG)	HFA 42 dB
-------------------------------------	-----------

Reference test acousto-mechanical sensitivity	
---	--

level (re. 1 μ N / 20 μ Pa)	HFA 35 dB
-------------------------------------	-----------

(60 dB SPL input, RTS)*	
-------------------------	--

Frequency range*	200–9850 Hz
------------------	-------------

Equivalent input noise*	24 dB SPL
-------------------------	-----------

Battery current*	1.8 mA
------------------	--------

(65 dB SPL input at 1 kHz)	
----------------------------	--

Total Harmonic Distortion*

70 dB SPL at 500 Hz	< 3%
---------------------	------

70 dB SPL at 800 Hz	< 0.3%
---------------------	--------

65 dB SPL at 1600 Hz	< 0.3%
----------------------	--------

60 dB SPL at 3200 Hz	< 0.3%
----------------------	--------

* Performed with reference test setting (RTS) of the gain control according to IEC 60118-9:2019. The RTS is set to full-on-gain minus 7 dB.

FOG = Full On Gain

RTS = Reference Test Setting

HFA = High Frequency Average

SPL = Sound Pressure Level re. 20 μ Pa

Environmental conditions

Condition	Minimum	Maximum
Operating temperature	+5 °C (41 °F)	+40 °C (104 °F)
Operating humidity	10% RH	90% RH
Operating pressure	700 hPa	1060 hPa
Transport temperature*	-10 °C (14 °F)	+55 °C (131 °F)
Transport humidity*	20% RH	95% RH
Storage temperature	+15 °C (59 °F)	+30 °C (86 °F)
Storage humidity	20% RH	90% RH

* Transport conditions require transport packaging used for sound processor.

Note

Battery performance deteriorates in temperatures below +5 °C.

Environmental protection

Your sound processor contains electronic components subject to the Directive 2012/19/EU on waste electrical and electronic equipment.

Help protect the environment by not disposing of your sound processor or batteries with your unsorted household waste. Please recycle your device, batteries and electronic items in according to your local regulations.

Magnetic Resonance Imaging (MRI)

The sound processor and other external accessories should never be brought into a room with an MRI machine, as damage to the sound processor or the MRI equipment could occur. The sound processor must be removed before entering a room where an MRI scanner is located.

If you are to undergo an MRI procedure, refer to the MRI Reference Card included in the document pack.

Electromagnetic compatibility (EMC)

Interference may occur in the vicinity of equipment marked with the following symbol:

Devices such as airport metal detectors, commercial theft detection systems, and Radio Frequency ID (RFID) scanners may produce strong electromagnetic fields. Some Baha users may experience a distorted sound sensation when passing through or near one of these devices. If this occurs, you should turn off the sound processor when in the vicinity of one of these devices. The materials used in the sound processor may activate metal detection systems. For this reason, you should carry the Security Control MRI Information Card with you at all times.

 Warning

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 in.) to any part of your sound processor, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

 Warning

Use of accessories, transducers and cables other than those specified or provided by Cochlear could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Regulatory information

Not all products are available in all markets. Product availability is subject to regulatory approval in the respective markets.

Equipment classification and compliance

Your sound processor is internally powered equipment Type B applied part as described in the international standard IEC 60601- 1:2005/A1:2012, Medical Electrical Equipment– Part 1: General Requirements for Basic Safety and Essential Performance.

This device complies with part 15 of the FCC (Federal Communications Commission) Rules. This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions:

- This device may not cause interference.
- This device must accept any interference, including interference that may cause undesired operation of the device.

Changes or modifications made to this equipment not expressly approved by Cochlear Bone Anchored Solutions AB may void the FCC authorization to operate this equipment. This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the recipient is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet or a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

FCC ID: QZ3BAHA6MAX

IC: 8039C-BAHA6MAX

HVIN: Baha 7

FVIN: 7.1.0

PMN: Cochlear Baha 7 Sound Processor

The model is a portable radio transmitter and receiver. It is designed not to exceed the emission limits for exposure to radio frequency (RF) energy set by the FCC and ISED. The safety distance of use is 0 mm.

The sound processor is designed not to exceed the emission limits according to CAN ICES-003 (B)/NMB-003(B).

Certification and applied standards

The products are in compliance with the following regulatory requirements:

- In EU: the device conforms to the General Safety and Performance Requirements according to Annex I of Medical Device Regulation 2017/745 (MDR) and essential requirements and other relevant provisions of Directive 2014/53/EU (RED). The declaration of conformity may be consulted at www.cochlear.com.
- Other identified applicable international regulatory requirements in countries outside the EU and US. Please refer to local country requirements for these areas.

Key to symbols

The following symbols may appear on your processor, accompanying accessories and/or packaging:

Consult instruction for use

Audio signal

Manufacturer

Keep dry

Catalogue number

CE mark and Notified
Body number

Serial number

UK Conformity mark with
approved body number

Unique device identifier

Medical device

Batch code

Rx Only

By prescription

Date of manufacture

Recyclable packaging

Temperature limitation

Waste of Electrical and Electronic Equipment

Note

Important information or advice.

Type B applied part

MR Unsafe

Interference risk

Caution (no harm)

Special care to be taken to ensure functionality. Could cause damage to equipment.

MR Conditional

Radio certification symbol for Korea

Warning (harmful)

Potential safety hazards and serious adverse reactions. Could cause harm to person.

ACMA symbol (Australian Communications and Media Authority)

IP42

IP 42 Ingress Protection Rating, protected against

- ingress of solid foreign objects
- dripping water

Radio certification symbol for Brazil

Radio certification symbol for Japan

Refer to instructions/ booklet

Trademark legal notice

For information regarding the compatibility of Cochlear's Sound Processors with Apple or Android devices, visit www.cochlear.com/compatibility. Cochlear's range of Smart Apps and Support Apps are available on App Store and/or Google Play. For information regarding the sound processors, operating systems and devices that are compatible with the Cochlear's range of Smart Apps and Support Apps, visit www.cochlear.com/compatibility.

Cochlear, Baha, 科利耳, コクレア, 코클리어, Hear now. And always, SmartSound, the elliptical logo, and marks bearing an ® or ™ symbol, are either trademarks or registered trademarks of the Cochlear group of companies (unless otherwise noted).

Android, Google Play and the Google Play logo are trademarks of Google LLC.

Apple, the Apple logo, iPhone, iPad and iPod are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Use of the Made for Apple badge means that an accessory has been designed to connect specifically to the Apple product(s) identified in the badge and has been certified by the developer to meet Apple performance standards. Apple is not responsible for the operation of this device or its compliance with safety and regulatory standards.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Cochlear Limited is under license.

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[GB] UK Responsible Person: Cochlear Europe Ltd, 6 Dashwood Lang Road, Bourne Business Park, Addlestone, Surrey, KT15 2HU, United Kingdom. Tel:+44 1932 26 3400

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[SG] Cochlear Limited (Singapore Branch), 238A Thomson Road #25-06, Novena Square Office Tower A, Singapore 307684, Singapore. Phone: +65 65533814

[CN] Cochlear Medical Device (Beijing) Co Ltd, Unit 2608-2617, 26th Floor, No.9 Building, No.91 Jianguo Road, Chaoyang District, Beijing 100022, P.R. China. Tel:+86 10 5909 7800

[HK] Cochlear (HK) Ltd, Room 1404-1406, 14/F, Leighton Centre, 77 Leighton Road, Causeway Bay, Hong Kong. Tel:+852 2530 5773

[IN] Cochlear Medical Device Company India PVT Ltd, Platina Bldg, Ground Floor, Plot No. C 59, G Block, BKC, Bandra East, Mumbai 400051 India. Tel:+91 22 6112 1111

[CO] Cochlear Colombia, Avenida Carrera 9 #115-06 Of. 1201 Edificio Tierra Firme, Bogota D.C., Colombia. Tel:+57 315 339 7169 / +57 315 332 5483

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Android, Google Play e o logotipo do Google Play são marcas comerciais da Google LLC.
[SE] Cochlear Bone Anchored Solutions AB, Konstruktionsvägen 14, 435 33 Mölnlycke, Sweden. Tel:+46 31 792 44 00
[CHREP] [CH] Cochlear AG, Peter Merian-Weg 4, 4052 Basel, Switzerland. Tel:+41 61 205 8204
[DE] Cochlear Deutschland GmbH & Co. KG, Mailänder Straße 4a, 30539 Hannover, Germany. Tel:+49 511 542 770
[SE] Cochlear Nordic AB, Konstruktionsvägen 14, 435 33 Mölnlycke, Sweden. Tel:+46 31 335 14 61
[CA] Cochlear Canada Inc, 2500-120 Adelaide Street West, Toronto, ON M5H 1T1, Canada. Tel:+1 800 483 3123