



Cochlear™ Osia® OSI200 Implant template

INSTRUCTIONS FOR USE

EN

FOR PROFESSIONALS

Contents

Intended purpose	4
Medical Indications	4
Contraindications	4
Intended patient population	4
Clinical benefits	4
Intended users	4
Serious incidents	5
Using the template	5
Storage	5
Disposal	5
Summary of safety and clinical performance	5
Labelling symbols	6
Warnings	7

Cochlear™ Osia®

OSI200 Implant template

Intended purpose

The OSI200 Implant template is intended to be used throughout the surgery to correctly position and attach the OSI200 Implant or OSI300 Implant.

Medical Indications

The OSI200 Implant template has no medical indication in itself. It will support the medical indications of compatible OSI Implant (OSI200 Implant and OSI300 Implant).

Contraindications

The OSI200 Implant template is contraindicated for use with Implants other than the OSI200 Implant or the OSI300 Implant.

Intended patient population

Refer to the OSI200 Implant Physician's Guide or the OSI300 Implant Physician's Guide.

Clinical benefits

The OSI200 Implant template functions as part of the Cochlear Osia System. Most recipients of a bone conduction hearing solution will experience an improved hearing performance and quality of life compared to unaided listening.

Intended users

The intended users of the OSI200 Implant template in a surgical environment are surgeons and qualified medical professionals (e.g. registered nurses).

Serious incidents

Serious incidents are rare, any serious incident in relation to your device should be reported to your Cochlear representative and to the medical device authority in your country, if available.

Using the template

The sterile single use OSI200 Implant template can be used:

- Before making the incision to mark the implant shape
- After making the incision, to check the implant position
- During the coil pocket creation to check for its size

Note: When using the OSI200 Implant template in a non-sterile area, make sure to use a new sterile implant template in the sterile field.

For details on the above procedures refer to the OSI200 Implant Physician's Guide or OSI300 Implant Physician's Guide.

Storage

For long term storage, store at room temperature. Keep dry. For details on storage conditions or technical specifications, refer to the OSI200 Implant Physician's Guide or OSI300 Implant Physician's Guide.










Disposal

OSI200 Implant templates that have been in patient contact should be placed into the correct clinical waste container for disposal. Follow the legal provisions for your country and the hygiene instructions of your hospital or clinic.

Summary of safety and clinical performance

<https://ec.europa.eu/tools/eudamed>

Labelling symbols

By prescription		Rx Only	
Date of manufacture		Manufacturer	
Authorised representative in the European Community		Use by date	
Warning. Potential safety hazards and serious adverse reactions. Could cause harm to person.			
UK Conformity mark with approved body number		Keep dry	
CE registration mark with notified body number		Medical device	

Warnings

For temporary use only. Not for implantation.

Supplied sterile. Sterilized using irradiation.

STERILE R

Sterile barrier



Single sterile barrier system with protective packaging outside



Do not reuse.
Do not resterilise.



Do not use if package is damaged and consult instructions for use



Do not use if template becomes non-sterile e.g. if dropped or mishandled in theatre after removal from packaging.

Compatible with OSI200 Implant and OSI300 Implant.

Not all products are approved in all markets.

AU Cochlear Ltd (ABN 96 002 618 073)
1 University Avenue, Macquarie University, NSW 2109,
Australia
Tel: +61 2 9428 6555

ECREP DE Cochlear Deutschland GmbH & Co. KG
Mailänder Straße 4 a, 30539 Hannover, Germany
Tel: +49 511 542 770

CHREP CH Cochlear AG
Peter Merian-Weg 4, 4052 Basel, Switzerland
Tel: +41 61 205 8204

US Cochlear Americas
10350 Park Meadows Drive, Lone Tree, CO 80124, USA
Tel: +1 (800) 523 5798

CA Cochlear Canada Inc
2500-120 Adelaide Street West, Toronto, ON M5H 1T1, Canada
Tel: +1 (800) 523 5798

GB UK Responsible Person: Cochlear Europe Ltd
6 Dashwood Lang Road, Bourne Business Park, Addlestone,
Surrey KT15 2HJ, United Kingdom
Tel: +44 1932 26 3400

BE Cochlear Benelux NV
Schaliënhoedreef 20 i, B-2800 Mechelen, Belgium
Tel: +32 15 79 55 11

FR Cochlear France S.A.S.
135 Route de Saint-Simon, 31035 Toulouse, France
Tel: +33 5 34 63 85 85 (International) or 0805 200 016
(National)

IT Cochlear Italia S.r.l.
Via Trattati Comunitari Europei 1957-2007 n.17,
40127 Bologna (BO), Italy
Tel: +39 051 601 53 11

SE Cochlear Nordic AB
Konstruktionsvägen 14, 435 33 Mölnlycke, Sweden
Tel +46 31 335 14 61

TR Cochlear Tıbbi Cihazlar ve Sağlık Hizmetleri Ltd. Şti.
Küçükbakkalköy Mah, Defne Sok, Büyükhanlı Plaza No:3 Kat:3
Daire: 9-10-11-12, 34750, Ataşehir, İstanbul, Türkiye
Tel: +90 216 538 5900

HK Cochlear (HK) Limited
Room 1404-1406, 14/F, Leighton Centre, 77 Leighton Road,
Causeway Bay, Hong Kong
Tel: +852 2530 5773

KR Cochlear Korea Ltd
2nd Floor, Yongsan Centreville Asterium, 25,
Hangang-daero 30 gil, Yongsan-gu, Seoul, Korea (04386)
Tel: +82 2 533 4450

CN Cochlear Medical Device (Beijing) Co., Ltd
Unit 2608-2617, 26th Floor, No.9 Building, No.91 Jianguo
Road,
Chaoyang District, Beijing 100022, P.R. China
Tel: +86 10 5909 7800

IN Cochlear Medical Device Company India Pvt. Ltd.
Ground Floor, Platina Building, Plot No C-59, G-Block,
Bandra Kurla Complex, Bandra (E), Mumbai – 400 051, India
Tel: +91 22 6112 1111

JP 株式会社日本コクレア(Nihon Cochlear Co Ltd)
〒113-0033 東京都文京区本郷2-3-7 お茶の水元町ビル
Tel: +81 3 3817 0241

AE Cochlear Middle East FZ-LLC
Dubai Healthcare City, Al Razi Building 64, Block A, Ground
Floor, Offices IR1 and IR2, Dubai, United Arab Emirates
Tel: +971 4 818 4400

PA Cochlear Latinoamérica S.A.
International Business Park, Building 3835, Office 403,
Panama Pacifico, Panama
Tel: +507 830 6220

NZ Cochlear NZ Limited
Level 4, Takapuna Towers, 19-21 Como St, Takapuna,
Auckland 0622, New Zealand
Tel: + 64 9 914 1983

www.cochlear.com



Cochlear, Osia, 科利耳, コクレア, 코클리어, Hear now. And always, SmartSound, the elliptical logo, and marks bearing an ® or ™ symbol, are either trademarks or registered trademarks of the Cochlear group of companies (unless otherwise noted).

© Cochlear Limited 2025.



P2323085 D2328711-V1 2025-03

D2328711-V1