Cochlear™ Baha®

SP Magnets

INSTRUCTIONS FOR USE

EN-GB	FI	KO	RU
AR	FR	LT	SK
BG	HE	LV	SL
CS	HR	MS	SR
DA	HU	NL	SV
DE	ID	NO	TH
EL	IS	PL	TR
ES	IT	PT	UK
ET	JA	PT-BR	ZH-CN
FA	KK	RO	ZH-TW

Contents

EN-GB.....	3	KO.....	86
AR.....	7	LT.....	90
BG.....	11	LV.....	94
CS.....	15	MS.....	98
DA.....	19	NL.....	102
DE.....	23	NO.....	107
EL.....	28	PL.....	111
ES.....	33	PT.....	115
ET.....	37	PT-BR.....	119
FA.....	41	RO.....	124
FI.....	46	RU.....	128
FR.....	50	SK.....	133
HE.....	54	SL.....	137
HR.....	58	SR.....	141
HU.....	62	SV.....	145
ID.....	66	TH.....	149
IS.....	70	TR.....	153
IT.....	74	UK.....	157
JA.....	78	ZN-CN.....	161
KK.....	82	ZH-TW.....	164

EN-GB Cochlear™ Baha®

SP Magnets Instructions for use

Intended use: The Cochlear Baha System uses bone conduction to transmit sounds to the cochlea (inner ear) with the purpose of enhancing hearing. The SP Magnet is intended to be used as part of the Cochlear Baha System to function as retention and coupling mechanism between the Baha implant magnet and the Baha Sound Processor.

Indications: The Cochlear Baha System is indicated for patients with conductive, mixed hearing loss and single-sided sensorineural deafness (SSD).

Clinical benefit: Most recipients of a bone conduction hearing solution will experience improved hearing performance and quality of life compared to unaided listening.

Support literature and hands-on workshops are available from Cochlear. Go to www.cochlear.com for contact information.

Cochlear accepts no responsibility for any adverse outcomes if used with products not recommended by Cochlear. Close cooperation in an interdisciplinary team is essential for a successful outcome.

Patients should have sufficient bone quality and quantity to support successful implant placement.

The SP Magnets are available in different strengths. If the magnet is too weak, the sound processor may fall off and if it is too strong, the patient may feel discomfort or experience skin

soreness and irritation. If the SP Magnet causes discomfort, first try a weaker magnet, if that does not help, attach one extra SP Magnet softpad.

Soft tissue thickness	Start SP Magnet strength
3 mm	1
4 mm	2
5 mm	3
6 mm	4

Contraindications:

- a) A pre-operative evaluation determines possible risks from the implant procedure and healing capacities of the bone and/or associated soft tissue.
- b) Insufficient bone quality or quantity.
- c) Patients with insufficient soft tissue thickness.

Children can be fitted, but treatment at an early age is dependent on the surgeon's evaluation regarding adequate bone depth and soft tissue thickness. In the United States and Canada, the placement of a bone anchored implant is contraindicated in children below the age of 5.

Possible complications: A 100 % rate of success cannot be guaranteed. Potential causes for failure of osseointegration include lack of adequate bone quantity/quality, trauma, infection, generalised diseases and surgical complications. Patients should be educated to contact a physician in case of significant trauma,

infection, skin irritation, pain, discomfort or implant loss. If the skin has been reduced and grows back to its original thickness, skin reduction surgery may be indicated.

Precautions: Do not fit an SP Magnet if the surgical incision is not sufficiently healed.

The sound processor and SP Magnet must be removed before entering a room where an MRI scanner is located and for any treatment involving ionizing radiation, such as for example X-ray, CT-scan and radiation therapy.

The implanted magnet may affect diagnostic information/ treatment effectiveness in the area around the magnet during treatment with ionizing radiation (e.g. X-Ray, CT-scan, radiation therapy). Take this into consideration before any such treatment. Counsel the patient or patient's caregivers that the Baha Attract System contains small parts that may be hazardous if swallowed or may cause choking if ingested or inhaled.

The Baha Attract System contains magnets that should be kept away from life supporting devices, e.g. cardiac pacemakers and ICDs (implantable cardioverter defibrillators) and magnetic ventricular shunts, as the magnets may affect the function of these devices. Contact the manufacturer of the specific device to find out more. Only Cochlear products should be used with each other. Combining products that are not dimensioned for mating can lead to mechanical failure, damage to tissue or to compromised audiological results.

Contains small parts that may present a choking hazard or suffocation.

Cleaning: Remove the soft pad and clean the SP Magnet with an alcohol wipe (70% ethanol) before use on new patient.

Warranty: We, Cochlear, warrant to you, the consumer of the Product, that: a. each Product is of merchantable quality; b. each Product is reasonably fit for the purpose or purposes for which it is supplied by us; and c. each Product will be free from defects in design, workmanship and materials for the Warranty Period. This warranty excludes liability for the defects or damage arising from, associated with, or related to the use of this Product with any non-Cochlear processing unit and/or any non-Cochlear implant. Warranty period for SP magnets: 2 years starting from the earliest date of either: the first fit date OR 3 months after the date shipped by Cochlear.

Serious incidents: Serious incidents are rare, any serious incident in relation to your device should be reported to your Cochlear representative and to the medical device authority in your country, if available.

 MD	Medical device	 UDI	Unique device identifier
	Keep dry	Store at room temperature	Store at room temperature

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**GB UK Responsible Person:
Cochlear Europe Ltd**,
6 Dashwood Lang Road,
Bourne Business Park, Addlestone,
Surrey, KT15 2HJ, United Kingdom.
Tel:+44 1932 26 3400

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www.cochlear.com

This material is intended for health professionals. If you are a consumer, please seek advice from your health professional about treatments for hearing loss. Outcomes may vary, and your health professional will advise you about the factors which could affect your outcome. Always read the instructions for use. Not all products are available in all countries. Please contact your local Cochlear representative for product information.

Contains small parts
that may present a
choking hazard.

Medical
device

Rx Only

Cochlear, Baha, 科利耳, コクレア, 코클리어, Hear now. And always, SmartSound, the elliptical logo, and marks bearing an ® or ™ symbol, are either trademarks or registered trademarks of the Cochlear group of companies (unless otherwise noted).

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CHREP CH Cochlear AG,
Hizmetleri Ltd. Sti.**, Küçükbakkalköy Mah,