



Cochlear[®]
Hear now. And always

Quick guide for Android™ connectivity

Cochlear™ Baha® 6 Max Sound Processor

Benefits of direct-to-device connectivity from compatible Android™ smartphones (ASHA).

- Stream your favorite movies, TV shows, music and more directly to your sound processor
- Stream phone calls directly to the sound processor
- Use the Cochlear™ Baha® Smart App

Compatibility

For up-to-date compatibility, please refer to www.cochlear.com/compatibility



android 

Pairing and connecting the sound processor

Below are example screens that shows how to pair your Baha 6 Max Sound Processor with an Android smartphone.

Please refer to the user guide for your compatible Android smartphone and follow the instructions on how to pair a Bluetooth® device.

1. Go to **Settings > Connections > Bluetooth** and make sure Bluetooth is turned on (fig. 1,2).
2. Restart the sound processor by opening and closing the battery door (fig. 3).




fig. 1




fig. 2



fig. 3



Cochlear®
Hear now. And always

Quick guide for Android connectivity

Cochlear™ Baha® 6 Max Sound Processor

Pairing and connecting the sound processor

3. Your sound processor will be listed under Available devices (fig. 4). If your sound processor does not appear, restart it and also tap 'Scan' in the top-right corner of your smartphone screen (fig. 5).
4. Tap on the name of your sound processor in the list to establish a connection (fig. 6).
5. The smartphone will pair and connect to the sound processor (fig. 7).




fig. 4




fig. 5

Unpairing your sound processor

To remove a previously paired sound processor from your smartphone, go to **Settings > Bluetooth** and tap on the cog wheel next to your sound processor (fig 7).

In the bottom-right corner of the screen that opens up (fig. 8), tap to 'Unpair'.

Your device will still show in the list of available Bluetooth devices.

To pair to it again, simply tap it.

To clear the list of available devices, turn Bluetooth off, then on.




fig. 6




fig. 7




fig. 8



Cochlear[®]
Hear now. And always

Quick guide for Android connectivity

Cochlear™ Baha® 6 Max Sound Processor

Streaming phone calls

Once paired, phone calls can be streamed directly to the connected sound processor.

When the sound processor connects, it should become the default audio streaming destination.

When you make a phone call, a menu automatically pops up allowing you to select where audio should be streamed.

When you receive a phone call (or during a phone call), you can press the “Bluetooth” icon (fig. 9) on the screen and select where audio should be streamed.




fig. 9

Re-connecting to your smartphone

The device will disconnect from your sound processor when you turn off the sound processor, when the sound processor battery level is too low, or when you are out of range.

To automatically reconnect, make sure the sound processor has sufficient battery power and your sound processor is within range of the device. You should hear a series of beeps in the sound processor confirming the connection.

If the connection is not automatically re-established, you may need to unpair your sound processor and repeat the pairing sequence.

Pairing to bilateral sound processors

To stream audio to two sound processors simultaneously they need to be bilaterally paired.

If you don't know if your sound processors are bilaterally paired, contact your and ensure that your two sound processors have been configured as a left and a right device and that Control Sync is ON.



Cochlear[®]
Hear now. And always

Quick guide for Android connectivity

Cochlear[™] Baha[®] 6 Max Sound Processor

Pairing to more than one smartphone

You can pair up to five different smartphones to your sound processor, but you can only be connected to one sound processor (or bilateral pair) at a time.

To manually change connection between devices, turn off Bluetooth on the connected device, and turn on Bluetooth on the device you want to connect to. It will then automatically reconnect.

Support

For more information and support on the Baha 6 Max Sound Processor and Android streaming, please contact your hearing care professional, refer to the sound processor user manual or visit the Cochlear website: www.cochlear.com

For further support visit www.cochlear.com/customer-service to find contact details to your local Cochlear Customer Service team.

Please seek advice from your health professional about treatments for hearing loss. Outcomes may vary, and your health professional will advise you about the factors which could affect your outcome. Always read the instructions for use. Not all products are available in all countries. Please contact your local Cochlear representative for product information.

Android is a trademark of Google LLC. The Android robot is reproduced or modified from work created and shared by Google and used according to terms described in the Creative Commons 3.0 Attribution License. Google Play and the Google Play logo are trademarks of Google LLC.

The Bluetooth[®] word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Cochlear Limited is under license.

©Cochlear Limited 2021. All rights reserved. Hear now. And always and other trademarks and registered trademarks are the property of Cochlear Limited or Cochlear Bone Anchored Solutions AB.