

Cochlear[™] Nucleus[®] Implants Magnetic Resonance Imaging (MRI) Guidelines

Asia Pacific

Hear now. And always



About this guide

This guide applies to Cochlear™ Nucleus® implants. It is intended for:

- Specialised health care professionals who prepare and perform MR scans
- Physicians who refer a Cochlear Nucleus implant recipient for an MR scan
- Cochlear Nucleus implant recipients and/or their carers

This guide provides information about the safe application of an MR scan on Cochlear Nucleus implant recipients.

MR scans performed under different conditions than those presented in this guide may result in severe patient injury or device malfunction.

Due to the risks associated with using MRI with an implanted medical device, it is important to read, understand, and comply with these instructions to prevent potential harm to the patient and/or device malfunction.

This guide should be read in conjunction with the relevant documents that accompany a Cochlear Nucleus implant, such as the Physician's Guide and Important Information Booklet. For more information, visit www.cochlear.com/warnings.

Symbols used in this guide



Note

Important information or advice.



Caution (no harm)

Special care to be taken to ensure safety and effectiveness.

Could cause damage to equipment.



Warning (harmful)

Potential safety hazards and serious adverse reactions.

Could cause harm to person.

Contents

- About this guide.....1
 - Symbols used in this guide 2
- MRI safety information..... 5
 - Bilateral recipients..... 5
 - Identifying the Cochlear Nucleus implant..... 5
 - X-ray information for identification of Cochlear Nucleus implants..... 6
 - X-ray guidelines 6
- Implant model identification 7
 - Cochlear Nucleus CI600 Series and CI500 Series implants..... 7
 - Cochlear Nucleus CI24RE Series, CI24R Series, CI24M Series and CI22M Series implants 9
- MRI safety information for Cochlear Nucleus implants 12
 - Implant magnet conditions for MRI..... 12
 - Indications for using MRI safely 14
 - CI600 Series implants 16
 - CI500 Series implants 17
 - CI24RE Series implants..... 18
 - CI24R and CI24M Series implants..... 19
 - CI22M Series implants..... 20
- Image interference and artefacts 21
- Preparation prior to an MRI examination 25
 - Cooperation between specialists 25
 - Considerations for implant magnet removal..... 26

Considerations for conducting an MRI examination.....	28
Prerequisites.....	28
Patient positioning	29
Patient comfort.....	29
Perform the MR scan	30
Performing an MR scan on other body locations.....	30
Cochlear Nucleus Implant Bandage and Splint Kit for MRI (MRI Kit) ..	31
Intended use	31
MRI Kit contraindications.....	32
Obtaining an MRI Kit	32
MRI Kit contents	32
Using the MRI Kit.....	33
Considerations after an MRI examination.....	38
With the implant magnet in place.....	38
With the implant magnet removed.....	38
Disposal	38
Considerations for referring physicians.....	39
Risks associated with MRI and Cochlear Nucleus implants	41
Labelling symbols	42
Legal statement	43
Trademark legal notice	43

MRI safety information

In order to determine if a patient may receive an MR scan, you must first identify the patient's Cochlear Nucleus implant model.

After you have identified the implant model, see *MRI safety information for Cochlear Nucleus implants* on page 12 to locate the MRI safety information for that specific implant model.



All external components of the Cochlear implant system (e.g. sound processors, remote assistants and related accessories) are MR Unsafe. The patient must remove all external components of their Cochlear implant system before entering a room where an MRI scanner is located.

Bilateral recipients

If one or more of the implants is a CI22M cochlear implant without a removable magnet, MRI is contraindicated.

If a bilateral recipient has implant models (other than the CI22M cochlear implant without a removable magnet), read the MRI safety information for each implant model relevant to the recipient. Use the MRI safety information of the recipient's implant model with the most restrictive MRI exposure requirements.

Identifying the Cochlear Nucleus implant

The implant model can be found on the patient's Cochlear patient identification card.

If the patient does not have their patient identification card with them, the implant type and model can be identified without surgical intervention. See *X-ray information for identification of Cochlear Nucleus implants* on page 6 and *Implant model identification* on page 7.

X-ray information for identification of Cochlear Nucleus implants

Cochlear Nucleus implants are made of metal and implanted under the skin behind the ear.




Figure 1: Location behind the ear for Cochlear Nucleus implants

X-ray guidelines

Lateral X-ray at 70 kV/ 3 mAs provides sufficient contrast to identify the implant.

A modified Stenver's view is not recommended for implant identification as implants may appear oblique.

Imaging should include an unobstructed view of antenna coils and implant bodies.

Bilateral recipients may have different implant models on either side of the head. A lateral skull X-ray with a 15 degree cranial tube angle will offset the implants in the image, enabling identifying features to be distinguished.

Implant model identification

Identifying features on Cochlear Nucleus implant X-ray images are explained in the following pages. Other implant models may have other identifying features.

Cochlear Nucleus CI600 Series and CI500 Series implants

Cochlear Nucleus CI600 Series implants - CI612, CI622 and CI632 and CI500 Series implants – CI512, CI522 and CI532 – do not have radiopaque characters.

Using an X-ray, CI500 Series and CI600 Series implants can be identified by the implant shape and electronic assembly layout. If further implant details are required, contact your Cochlear representative who will provide instructions on how to determine the following:

- Manufacturer
- Model
- Year of manufacture

*Not all products are available in all countries.
Please contact your local Cochlear representative for product information.

The electronic assembly layout is identical for Cochlear CI600 and CI500 Series implants. The unique identifier for CI600 Series implants is the magnet shape and the three holes next to the magnet, as illustrated in the table below.



CI600 Series implant X-ray	CI500 Series implant X-ray	Unique identifier
		<ol style="list-style-type: none"> 1. Three holes adjacent to magnet 2. Magnet shape 3. Round shape at coil exit end of electronic assembly layout 4. Series of wire connectors that are visible on both sides of the electronic assembly 5. Square implant body shape

Table 1: CI600 & CI500 Series implant identified by their shape and electronic assembly

Cochlear Nucleus CI24RE Series, CI24R Series, CI24M Series and CI22M Series implants

Cochlear Nucleus implants that can be identified by the radiopaque characters printed on them are:

- CI24RE Series: CI422, CI24REH, CI24RE (CA), CI24RE (CS) and CI24RE (ST)
- CI24R Series: CI24R (CA), CI24R (CS), CI24R (ST)
- CI24M Series: CI24M, CI 11+11+2M and ABI24M
- CI22M Series: CI22M

There are three sets of radiopaque characters printed on each implant.

1. The first character identifies the manufacturer – ‘C’ indicates Cochlear Ltd.
2. The second (middle) character identifies the implant model.
3. The third character indicates the year of manufacture. To determine the year of manufacture of your implant, contact your Cochlear representative.

Implant model	Location of second (middle) radiopaque character set	Radiopaque characters
CI422		13
CI24REH		6
CI24RE (CA)		5
CI24RE (CS)		7
CI24RE (ST)		4

Table 2: CI24RE Series implants identified by radiopaque characters

Implant model	Location of second (middle) radiopaque character set	Radiopaque characters
CI24R (CA)		2
CI24R (CS)		C
CI24R (ST)		H

Table 3: CI24R Series implants identified by radiopaque characters

Implant model	Location of second (middle) radiopaque character set	Radiopaque characters
CI24M		T
CI 11+11+2M		P
ABI24M		G

Table 4: CI24M Series implants identified by radiopaque characters

Implant model	Location of second (middle) radiopaque character set	Radiopaque characters
CI22M with removable magnet		L or J
CI22M without removable magnet		Z

Table 5: CI22M Series implants identified by radiopaque characters

MRI safety information for Cochlear Nucleus implants



Non-clinical testing has demonstrated that Cochlear Nucleus implants are MR Conditional.

Implant magnet conditions for MRI

For some implant models and MRI field strengths, bandaging with an MRI Kit is required, or the implant magnet needs to be surgically removed. Refer to the *Table 6* on page 13 for information on each Nucleus implant model.

Implant type	MRI field strength (T)	Remove implant magnet Yes/No	MRI Kit required Yes/No
CI600 Series implants			
CI612, CI622, CI632	1.5	No	No
	3		
CI500 Series implants			
CI512, CI522, CI532	1.5	No	Yes
	3	Yes	No
CI24RE Series implants			
CI422, CI24REH (Hybrid L24), CI24RE (CA), CI24RE (ST)	1.5	No	Yes
	3	Yes	No
CI24R and CI24M Series implants			
CI24R (CA), CI24R (CS), CI24R (ST), CI24M, ABI24M	1.5	No	Yes
	3	Yes	No
CI 11+11+2M	1.5	No	Yes
	3	MRI is contraindicated	
CI22M Series implants			
CI22M with removable magnet	1.5	No	Yes
	3	MRI is contraindicated	
CI22M without removable magnet	1.5	MRI is contraindicated	
	3		

Table 6: Implant magnet conditions for MRI

Indications for using MRI safely

The MRI safety information provided in these guidelines only applies to 1.5 T and 3 T MRI horizontal scanners (closed bore or wide bore) with a circularly polarised (CP) RF field for a maximum scan time of 60 minutes.



Warning

MR scans at 3 T must be performed in quadrature mode or CP mode for the radio frequency (RF) transmit coil. Using a multi-channel mode may result in localised heating above safe levels.

A patient with one or two of these devices can be safely scanned in an MR system meeting conditions on the following pages. All scans shall be performed according to the specified SAR limits for the relevant implant.

Consider the following prior to scanning:

- Determine if the magnet should be removed, or if the Cochlear Nucleus Implant Bandage and Splint Kit for MRI (MRI Kit) is required. See *Table 6: Implant magnet conditions for MRI* on page 13.
- Maximum allowable MRI scan time is 60 minutes of continuous scanning, with the SAR limitations provided in this guide. See *Indications for using MRI safely* on page 14.
- Transmit/receive head coils and whole body coils may be safely used within the recommended SAR limits. Refer to the MRI safety information and recommended SAR limit tables on the following pages.
- Local cylindrical transmit/receive coils may be safely used, without SAR restriction, provided that the distance between the entire implant and the end of the local RF coil is at least equal to the radius of the local RF coil.
- It is safe to use local cylindrical RF receive only coils with cochlear implants during MRI scanning, provided SAR limits for the transmit coil have not been exceeded.




Figure 2: Landmark locations

CI600 Series implants

CI600 Series implants can be safely scanned at least ten times without any adverse effect on magnet strength.

Implant type	MRI field strength (T)	Maximum allowable spatial gradient field (T/m)	Head average SAR limit (W/kg) Using transmit /receive head coil	Whole body average SAR limit (W/kg) Landmark location	
				<40 cm from top of head	≥40cm from top of head
CI612	1.5	20	<2	<1	<2
CI622					
CI632					
CI612	3	20	<1	<0.5	<1
CI622				<0.4	
CI632				<0.4	

Table 7: MRI safety information and recommended SAR levels for CI600 Series implants

CI500 Series implants

Implant type	MRI field strength (T)	Maximum allowable spatial gradient field (T/m)	Head average SAR limit (W/kg) Using transmit /receive head coil	Whole body average SAR limit (W/kg)	
				Landmark location	
				<40 cm from top of head	≥40cm from top of head
CI512	1.5	20	<2	<1	<2
CI522					
CI532					
CI512	3	20	<1	<0.5	<1
CI522				<0.4	
CI532				<0.4	

Table 8: MRI safety information and recommended SAR levels for CI500 Series implants

CI24RE Series implants

Implant type	MRI field strength (T)	Maximum allowable spatial gradient field (T/m)	Head average SAR limit (W/kg) Using transmit /receive head coil	Whole body average SAR limit (W/kg) Landmark location	
				<40 cm from top of head	≥40cm from top of head
CI422	1.5	20	<2	<1	<2
CI24REH					
CI24RE (CA)					
CI24RE (ST)					
CI422	3	20	<1	<0.5	<1
CI24REH					
CI24RE (CA)					
CI24RE ST)					

Table 9: MRI safety information and recommended SAR levels for CI24RE Series implants

CI24R and CI24M Series implants

Implant type	MRI field strength (T)	Maximum allowable spatial gradient field (T/m)	Head average SAR limit (W/kg) Using transmit /receive head coil	Whole body average SAR limit (W/kg) Landmark location	
				<40 cm from top of head	≥40cm from top of head
CI24R (CA)	1.5	20	<2	<1	<2
CI24R (CS)					
CI24R (ST)					
CI24M					
ABI24M					
CI 11+11+2M	1.5	20	<1	<0.5	<1
CI24R (CA)	3	20	<1	<0.5	<1
CI24R (CS)					
CI24R (ST)					
CI24M					
ABI24M					
CI 11+11+2M	3	MRI is contraindicated			

Table 10: MRI safety information and recommended SAR levels for CI24R and CI24M Series implants

CI22M Series implants

Implant type	MRI field strength (T)	Maximum allowable spatial gradient field (T/m)	Head average SAR limit (W/kg) Using transmit /receive head coil	Whole body average SAR limit (W/kg) Landmark location	
				<40 cm from top of head	≥40cm from top of head
CI22M with removable magnet	1.5	20	<2	<1	<2
	3	MRI is contraindicated			
CI22M without removable magnet	1.5	MRI is contraindicated			
	3				

Table 11: MRI safety information and recommended SAR levels for CI22M Series implants

Image interference and artefacts

The Cochlear Nucleus implant will create shadowing on the MR image near the implant, resulting in a loss of diagnostic information.

If inspecting near the implant, consider removing the implant magnet as MR image quality may be compromised with it in place.

If the implant magnet needs to be removed, refer the patient to an appropriate physician to arrange for the magnet to be removed before the MR scan.

The image artefact results¹ on the next page are based on worst-case spin echo scenarios. The images show maximum artefact extension from the centre of the implant and are representative of the axial results across all implants. The optimisation of scan parameters can be used to minimise the extent of the artefact. Tables detailing individual artefact sizes per implant model are on the following pages.

1. Image artefact testing was undertaken according to ASTM F2119 (Standard Test Method for Evaluation of MR Image Artefacts from Passive Implants) with worst case spin echo results provided.

For bilateral implant recipients, the image artefacts as shown below are mirrored on the opposite side of the head for each implant. There may be some extension of the artefact between the implants.


Implant magnet in place	Implant magnet removed
	
9.2 cm (3.6 in)	3.7 cm (1.5 in)

Table 12: Maximum artefact extension at 1.5 T across all implant types

Implant magnet in place (CI600 Series only)	Implant magnet removed
	
10.2 cm (4.0 in)	4.1 cm (1.6 in)

Table 13: Maximum artefact extension at 3 T across all implant types

	MRI field strength (T)	Maximum artefact radius [cm]	
		Implant magnet in place	Implant with non-magnetic cassette
		Axial	Axial
CI600 Series implants			
CI612, CI622, CI632	1.5	9.2	4.8
	3	10.2	5.6

Table 14: Artefact dimensions for CI600 Series implants

	MRI field strength (T)	Maximum artefact radius [cm]	
		Implant magnet in place	Implant magnet removed
		Axial	Axial
CI500 Series implants			
CI512, CI522, CI532	1.5	5.7	3.7
	3	N/A*	4.1
CI24RE Series implants			
CI422, CI24REH CI24RE (CA), CI24RE (ST)	1.5	5.5	3.2
	3	N/A*	3.4
CI24R Series implants			
CI24R (CA), CI24R (CS), CI24R (ST)	1.5	5.5	3.2
	3	N/A*	3.4
CI24M Series implants			
CI24M, ABI24M	1.5	5.5	3.7
	3	N/A*	4.1
CI 11+11+2M	1.5	5.5	3.7
	3	MRI is contraindicated	
CI22M Series implants			
CI22M with removable magnet	1.5	5.5	6
	3	MRI is contraindicated	
CI22M without removable magnet	1.5	MRI is contraindicated	
	3		

Table 15: Artefact dimensions for CI500, CI24RE, CI24R, CI24M and CI22M Series implants

* Surgically remove the implant magnet before MR scans at 3 T.

Preparation prior to an MRI examination

Cooperation between specialists

Preparing for and conducting an MRI examination for implant recipients requires cooperation between a specialist for the device and/or Cochlear Nucleus implant physician, referring physician and radiologist or MR technologist.

- **Cochlear Nucleus implant device specialist** – Knows the implant type and where to find the correct MR parameters for the implant.
- **Referring physician** – Knows the location of the MR scan and diagnostic information required, and makes a decision on whether the implant magnet needs to be removed for the MRI examination.
- **Cochlear Nucleus implant physician** – If requested by the referring physician, surgically removes the implant magnet and replaces with a non-magnetic plug or non-magnetic cassette. After the MR scan, the implant physician replaces it with a new sterile replacement implant magnet.
- **Radiologist and MR technologist** – Sets up the MR scan using the correct MR parameters and counsels the implant recipient during the MRI examination.

Considerations for implant magnet removal

If the implant magnet needs to be removed prior to an MRI examination, close coordination is required between the specialists to perform the implant magnet removal, MR scan, and subsequent implant magnet replacement.

For CI600 Series implant recipients, if single or multiple MRI examinations on the head are needed with the magnet removed, the implant magnet must be replaced (in a sterile surgical environment) with a non-magnetic cassette.



Warning

To prevent infection, do not leave the magnet pocket empty for CI600 implants. When removing the magnet cassette, replace the magnet cassette with a non-magnetic cassette.

For CI24RE, CI24R, CI24M, CI22M and CI500 Series implant recipients requiring multiple MRI examinations over a period of time, the implant magnet is removed and replaced with a sterile non-magnetic plug. If only a single MRI is required the magnet recess can remain empty. See *Implant magnet conditions for MRI* on page 12.

In the magnet's absence, the non-magnetic plug prevents fibrous tissue growing into the implant recess. Such growth would make implant magnet replacement difficult.



Caution

Non-magnetic plugs for CI500 Series implants are a different size to non-magnetic plugs for CI24RE Series implants. Ensure the correct plug is used.

With the non-magnetic cassette or non-magnetic plug in place, MR scans can be done at both 1.5 T and 3 T without the need for bandaging or use of the MRI Kit.



Note

While the magnet is removed, the recipient must wear a Cochlear Disk Retainer to hold their sound processor coil in place. Disk retainers are available from Cochlear.

When there is no further need for MRI examinations, the non-magnetic cassette or non-magnetic plug is removed and replaced by a new sterile replacement magnet.

The non-magnetic cassette, non-magnetic plug, replacement magnet cassette and sterile replacement magnet are supplied separately in sterile packs. All are single-use items.

Considerations for conducting an MRI examination

These guidelines are specific to Cochlear Nucleus implants and supplement other MRI examination considerations specified by the MRI machine manufacturer or protocols at the MRI facility.

Prerequisites

The following additional conditions must be met:

- The implant model has been identified. See *Implant model identification* on page 7.
 - For additional information on bilateral recipients, see *Bilateral recipients* on page 5.
- The artefact has been considered, and there is still diagnostic value in performing the MR scan. See *Image interference and artefacts* on page 21.
- If the referring physician has prescribed removal of the implant magnet for the MR scan, the implant magnet has been surgically removed. See *Preparation prior to an MRI examination* on page 25 and *Table 6: Implant magnet conditions for MRI* on page 13.
- The MRI Kit is required for MR scans at 1.5 T with the implant magnet in place for CI500 Series, CI24RE Series, CI24R Series, CI24M Series and CI22M Series implants.
 - See *Using the MRI Kit* on page 33 for instructions on how to apply the MRI Kit prior to the MR scan.
- Remove the sound processor before entering the MRI room. The sound processor is MR Unsafe.

Patient positioning

For safety, the patient should be in a supine position (lying flat on back, face upward) prior to entering the MRI bore.

Align the patient's head with the bore axis of the MRI machine. Advise the patient to lie as still as possible and to not move their head during the MR scan.

Caution

Ensure that the patient does not move more than 15 degrees (15°) from the centreline (Z-axis) of the bore during the MR scan.

Failure to position the patient correctly prior to the MR scan may result in increased torque on the implant and cause pain.

Patient comfort

Explain to the patient that they may sense the implant magnet moving. For implants where use of the MRI Kit is required, the MRI Kit will reduce the likelihood of the implant magnet moving. However, they may still sense resistance to movement as pressure on the skin. The sensation will be similar to pressing down firmly on the skin with the thumb.

If the patient experiences pain, consult the patient's physician to determine if the implant magnet should be removed or if a local anaesthetic may be applied to reduce discomfort.

Caution

If administering local anaesthetic, take care not to perforate the implant silicone.

In addition, explain to the patient that they may perceive sounds during the MR scan.

Perform the MR scan

The MR scan must be performed using the MRI safety information identified for the patient's implant model. See *Identifying the Cochlear Nucleus implant* on page 5 and the *MRI safety information and recommended SAR levels* tables in *MRI safety information for Cochlear Nucleus implants* on page 12.

Performing an MR scan on other body locations

When an implant recipient requires an MR scan on a location of their body away from the implant site, you must still follow the MRI safety information for the recipient's implant model. See *Identifying the Cochlear Nucleus implant* on page 5 and the *MRI safety information and recommended SAR levels* tables in *MRI safety information for Cochlear Nucleus implants* on page 12.

Cochlear Nucleus Implant Bandage and Splint Kit for MRI (MRI Kit)

Intended use

The Cochlear Nucleus Implant Bandage and Splint Kit for MRI (MRI Kit) is intended to be used on Cochlear Nucleus implant recipients to prevent implant magnet dislodgement during MR scans at 1.5 T.

The MRI Kit is for single-use only.

The MRI Kit is intended for use with the following Cochlear Nucleus implants for both unilateral and bilateral recipients:

- CI500 Series – CI512, CI522, and CI532
- CI24RE Series – CI422, CI24REH, CI24RE (CA), CI24RE (CS) and CI24RE (ST)
- CI24R Series – CI24R (CA), CI24R (CS) and CI24R (ST)
- CI24M Series – CI24M, CI 11+11+2M and ABI24M
- CI22M Series – CI22M (with removable magnet)

See *Table 6: Implant magnet conditions for MRI* on page 13 and *Bilateral recipients* on page 5 for full details on performing an MRI examination safely.



Warning

Although unlikely with the use of the MRI Kit, it is possible for the magnet to move during MRI and dislodge from the implant magnet pocket. In this case, surgical intervention to reposition or replace the magnet would be required.

MRI Kit contraindications

The MRI Kit is contraindicated for use with:

- CI22M Series – CI22M implants with non-removable magnet
- MR scans other than 1.5 T.

See the elasticised compression bandage labelling for related contraindications when using this product.

Obtaining an MRI Kit

Contact the nearest Cochlear office or official distributor to order an MRI Kit.

MRI Kit contents

Item	Description
Flat plastic splints	To be placed against the skin over the implant magnet site. For bilateral recipients, use one splint for each implant.
Elasticised compression bandage	For securing the splint against the implant magnet site.
Surgical tape	For securing the bandage and splint in place.

Using the MRI Kit

Follow this procedure to use the MRI Kit. When used as instructed, the supplied splint and bandage should reduce the likelihood of implant magnet movement when in or near the MRI scanner.

1. Preparation

1. Prior to entering the MRI room and before removing the sound processor, draw an outline of the BTE sound processor coil or OTE all-in-one unit on the patient's head - see *Figure 3* and *Figure 4*. Once the sound processor has been removed, mark the centre position of the outline; this is the implant magnet location. If necessary, shave the patient's head at the implant magnet location so this marking is more visible and easier to locate during the splinting process. This marking is essential to ensure that the splint is placed in the correct location. Repeat this step for bilateral recipients.



Note

Once the sound processor has been removed, the implant recipient may no longer be able to hear.




Figure 3: Location of the BTE (behind the ear) sound processor unit, sound processor coil and coil magnet




Figure 4: Location of the OTE (off the ear) sound processor unit, coil and external magnet

2. In the event that the location of the implant magnet has not been marked, it can be located by:
 - Using ferromagnetic material, such as a paper clip - the material will be attracted to the implant magnet.

⚠ Warning

The ferromagnetic material must be removed before entering the MRI room.

- Touch - gently feel around the implant site to locate the position of the implant coil. The implant is comprised of two components; the round implant coil and the implant body. See *Figure 5* below. The implant magnet will be at the centre of the implant coil.

2. Bandaging

1. Use a splint from the MRI Kit and centre it over the implant magnet site (as marked) against the skin. For bilateral recipients, use one splint for each implant. Ensure each splint is held in place over the implant magnet. See *Figure 5* below for the implant magnet location. You may need the assistance of another person to hold each splint in place while you bandage. Otherwise, use the supplied tape to maintain the splint position prior to bandaging.




Figure 5: Location of the implant magnet on CI500 Series (left side) and CI24RE Series (right side) implants

2. Use the elasticised compression bandage from the MRI Kit and ensure the centre line of the bandage is over the implant magnet site and the splint is fully covered. See **Figure 6** below.




Figure 6: Fitting the MRI Kit splint and compression bandage

3. Use a minimum of two bandage layers at full stretch (no elasticity remaining in the bandage). When the bandage is at its maximum tightness, the small rectangular tension markers will stretch to become square in shape. See *Figure 7* below.




Figure 7: Comparison of compression bandage tightness

4. Use the surgical tape from the MRI Kit to secure the bandage by wrapping two surgical tape layers around the head, over the bandage centre line. Ensure the tape ends overlap.
5. Conduct the MR scan.
6. Once the MR scan is complete, follow the instructions in *Considerations after an MRI examination* on page 38.

Considerations after an MRI examination

With the implant magnet in place

Remove the MRI Kit bandage and splint.

After the patient leaves the MRI room, ask the patient to place the sound processor on their head and turn it on.

Confirm:

- the placement of the sound processor is correct
- there is no discomfort
- sound is perceived as normal.

If there is discomfort or a change in sound perception, or problems with the placement of the sound processor, ask the patient to seek assistance from their implant clinician as soon as possible.

With the implant magnet removed

See *Considerations for implant magnet removal* on page 26.

Disposal

The MRI Kit can be disposed of as normal hospital or household waste, or in accordance with local regulations.

Considerations for referring physicians

If you are a physician referring a Cochlear Nucleus implant recipient for an MR scan, it is essential that you consider the following:

- Understand and inform the patient of the risks associated with MRI. See *Risks associated with MRI and Cochlear Nucleus implants* on page 41.
- Understand the conditions for an MR scan and ensure that there is a clear indication for the MRI examination. See *Implant magnet conditions for MRI* on page 12 and *Indications for using MRI safely* on page 14.
- Identify if the patient has any other medical device implants, active or abandoned. If another implanted device is present, verify MRI compatibility before conducting an MRI examination.

If MRI safety information for the implanted devices is not followed, the potential risks include:

- movement or damage to the device
 - weakening of the implant magnet
 - uncomfortable sensation for the patient
 - skin or tissue trauma for the patient.
- Cochlear has evaluated the interaction of implants described in this guide with other nearby implanted devices during MRI scanning, and there is no increased heating risk to the Cochlear implant.
 - The Cochlear Nucleus implant will create shadowing on the MR image in the vicinity of the implant, resulting in a loss of diagnostic information. Refer to the relevant artefact dimension tables in *Image interference and artefacts* on page 21.
 - For MR scans at 1.5 T or 3 T, identify if the implant magnet needs to be removed. See *Implant magnet conditions for MRI* on page 12.
 - For MR scans on a body location away from the implant site, MRI safety information for the recipient's implant model must still be followed. See *Performing an MR scan on other body locations* on page 30.



Figure 8: CI600 and CI500 Series implant with removable magnet

Consider the following:

- If the required diagnostic information is in the area of the implant, the implant magnet may need to be removed.
- Timing of the implant surgery and MRI exposure.
- Age and general health of the implant recipient and time to recover from the implant magnet surgery or potential trauma.
- Existing or potential for tissue scarring in the location of the implant magnet.
- If the implant magnet needs to be removed, refer the patient to an appropriate physician to arrange for the magnet to be removed before the MR scan.
- If the implant magnet is retained for an MR scan at 1.5 T, an MRI Kit must be obtained beforehand for use during the MR scan, except for CI600 Series implants. See *Obtaining an MRI Kit* on page 32.
- For CI600 Series cochlear implants, see *Table 6: Implant magnet conditions for MRI* on page 13 for full details on performing an MRI examination safely. For all other implants, please also see *Using the MRI Kit* on page 33.

Risks associated with MRI and Cochlear Nucleus implants

The potential risks of performing MRI examinations on patients with Cochlear Nucleus implants include:

- **Device movement**

Scanning outside of the parameters contained in these guidelines may lead to the implant magnet or device moving out of position during an MRI examination causing skin or tissue trauma.
- **Damage to the device**

MRI exposure beyond the values contained in these guidelines may cause damage to the device.
- **Weakening of implant magnet**
 - Scanning at static magnetic field strengths at values other than those contained in these guidelines may lead to a weakening of the implant magnet.
 - Incorrect patient positioning prior to the MR scan or head movement during the scan may result in implant magnet demagnetisation.
- **Uncomfortable sensation**

MRI exposure beyond the values contained in these guidelines may result in the patient perceiving sound or noise and/or pain.
- **Implant heating**

Use the recommended SAR values contained in these guidelines to ensure the implant does not heat beyond safe levels.
- **Image artefact**

The Cochlear Nucleus implant will create shadowing on the MR image in the vicinity of the implant, resulting in a loss of diagnostic information.

If inspecting near the implant, removal of the implant magnet should be considered as MR image quality may be compromised with it in place.

Labelling symbols

The following symbols may appear on the product, the components and/or the packaging.



Consult instructions for use



Refer to instruction manual



Specific warnings or precautions associated with the device, which are not otherwise found on the label



Manufacturer



Catalogue number



Authorised representative in the European Community



Keep dry



Do not use if package is damaged



Recyclable packaging

Rx Only

By prescription



MR Conditional

Legal statement

Cochlear implant systems are protected by one or more international patents. The statements made in this guide are believed to be true and correct as of the date of publication. However, specifications are subject to change without notice.

© Cochlear Limited 2021. All rights reserved.

Trademark legal notice

Cochlear implant systems are protected by one or more international patents.

ACE, Advance Off-Stylet, AOS, AutoNRT, Autosensitivity, Beam, Bring Back the Beat, Button, Carina, Cochlear, 科利耳, コクレア, 코클리어, Cochlear SoftWear, Codacs, Contour, コントウア, Contour Advance, Custom Sound, ESPrit, Freedom, Hear now. And always, Hugfit, Hybrid, Invisible Hearing, Kanso, MET, MicroDrive, MP3000, myCochlear, mySmartSound, NRT, Nucleus, Osia, Outcome Focused Fitting, Off-Stylet, Profile, Slimline, SmartSound, Softip, SPrint, True Wireless, the elliptical logo, and Whisper are either trademarks or registered trademarks of Cochlear Limited. Ardiium, Baha, Baha SoftWear, BCDrive, DermaLock, EveryWear, Human Design, Piezo Power, SoundArc, Vistafix, and WindShield are either trademarks or registered trademarks of Cochlear Bone Anchored Solutions AB.

Hear now. And always

Cochlear Ltd (ABN 96 002 618 073) 1 University Avenue, Macquarie University, NSW 2109, Australia
Tel: +61 2 9428 6555 Fax: +61 2 9428 6352

Cochlear Ltd (ABN 96 002 618 073) 14 Mars Road, Lane Cove, NSW 2066, Australia
Tel: +61 2 9428 6555 Fax: +61 2 9428 6352

ECOREP Cochlear Deutschland GmbH & Co. KG Karl-Wiechert-Allee 76A, 30625 Hannover, Germany
Tel: +49 511 542 770 Fax: +49 511 542 7770

Cochlear Americas 10350 Park Meadows Drive, Lone Tree, CO 80124, USA
Tel: +1 303 790 9010 Fax: +1 303 792 9025

Cochlear Canada Inc 2500-120 Adelaide Street West, Toronto, ON M5H 1T1, Canada
Tel: +1 416 972 5082 Fax: +1 416 972 5083

Cochlear AG EMEA Headquarters, Peter Merian-Weg 4, 4052 Basel, Switzerland
Tel: +41 61 205 8204 Fax: +41 61 205 8205

Cochlear Europe Ltd 6 Dashwood Lang Road, Bourne Business Park, Addlestone, Surrey KT15 2HJ, United Kingdom
Tel: +44 1932 26 3400 Fax: +44 1932 26 3426

Cochlear Benelux NV Schaliënhoedreef 20 i, B-2800 Mechelen, Belgium
Tel: +32 15 79 55 11 Fax: +32 15 79 55 70

Cochlear France S.A.S. 135 Route de Saint-Simon, 31035 Toulouse, France
Tel: +33 5 34 63 85 85 (International) or 0805 200 016 (National) Fax: +33 5 34 63 85 80

Cochlear Italia S.r.l. Via Trattati Comunitari Europei 1957-2007 n.17, 40127 Bologna (BO), Italy
Tel: +39 051 601 53 11 Fax: +39 051 39 20 62

Cochlear Nordic AB Konstruktionsvägen 14, 435 33 Mölnlycke, Sweden
Tel: +46 31 335 14 61 Fax: +46 31 335 14 60

Cochlear Tibbi Cihazlar ve Sağlık Hizmetleri Ltd. Şti.
Çubuklu Mah. Boğaziçi Cad., Boğaziçi Plaza No: 6/1, Kavacık, TR-34805 Beykoz-Istanbul, Turkey
Tel: +90 216 538 5900 Fax: +90 216 538 5919

Cochlear (HK) Limited Room 1404-1406, 14/F, Leighton Centre, 77 Leighton Road, Causeway Bay, Hong Kong
Tel: +852 2530 5773 Fax: +852 2530 5183

Cochlear Korea Ltd 1st floor, Cheongwon Building 33, Teheran-ro 8 gil, Gangnam-gu, Seoul, Korea
Tel: +82 2 533 4450 Fax: +82 2 533 8408

Cochlear Medical Device (Beijing) Co., Ltd
Unit 2608-2617, 26th Floor, No.9 Building, No.91 Jianguo Road, Chaoyang District, Beijing 100022, P.R. China
Tel: +86 10 5909 7800 Fax: +86 10 5909 7900

Cochlear Medical Device Company India Pvt. Ltd.
Ground Floor, Platina Building, Plot No C-59, G-Block, Bandra Kurla Complex, Bandra (E), Mumbai – 400 051, India
Tel: +91 22 6112 1111 Fax: +91 22 6112 1100

株式会社日本コクレア (Nihon Cochlear Co Ltd) 〒113-0033 東京都文京区本郷2-3-7 お茶の水元町ビル
Tel: +81 3 3817 0241 Fax: +81 3 3817 0245

Cochlear Middle East FZ-LLC
Dubai Healthcare City, Al Razi Building 64, Block A, Ground Floor, Offices IR1 and IR2, Dubai, United Arab Emirates
Tel: +971 4 818 4400 Fax: +971 4 361 8925

Cochlear Latinoamérica S.A.
International Business Park, Building 3835, Office 403, Panama Pacifico, Panama
Tel: +507 830 6220 Fax: +507 830 6218

Cochlear NZ Limited
Level 4, Takapuna Towers, 19-21 Como St, Takapuna, Auckland 0622, New Zealand
Tel: + 64 9 914 1983 Fax: 0800 886 036

www.cochlear.com