


Cochlear™ Baha® SoundArc



Component	Material	Additive/Treatment
Internal bow	Stainless steel spring wire, SUS 304	Electro Deposition Coating Powercron 6000CX, black
Connector disc	Aluminum AL6061-T6	Anodised black
Silicone Tip with pegs - black	MOMENTIVE Silopren LSR 4040	Pigment black CAS No. 1317-61-9 Poly(dimethylsiloxane) CAS No. 63148-62-9
Cover Tube – black	GRILAMID L 25 NATURAL 6112	PREBLACK PE 700-30

KEY DIMENSIONS

Sizes	X-Small	Small	Medium	Large	X-Large
Dimensions	30.6 cm	31.5 cm	33 cm	36 cm	39 cm

www.Cochlear.com/US

Cochlear Americas
13059 East Peakview Avenue
Centennial, CO 80111 USA
Telephone: 1 303 790 9010
Support: 1 800 483 3123

Cochlear Canada Inc.
2500-120 Adelaide Street West
Toronto, ON M5H 1T1 Canada
Support: 1 800 483 3123

©Cochlear Limited 2018. All rights reserved. Hear now. And always and other trademarks and registered trademarks are the property of Cochlear Limited or Cochlear Bone Anchored Solutions AB. The names of actual companies and products mentioned herein may be the trademarks of their respective owners.