

Cochlear™ Nucleus® Kanso® 2 Sound Processor Upgrade Order Form

PLEASE FILL OUT ALL FIELDS UNLESS INDICATED OTHERWISE.

01 Recipient Information

Name: _____

Email: _____

Phone: _____

Address: _____

Date of Birth: _____

Guardian Name (if child): _____

Gender: Male Female Non-binary Prefer not to say

Implant Side: Left Right

Preferred Language: English French

Additional comments (if needed): _____

02 Clinic Information

Clinic Name: _____

Clinic Email: _____

Audiologist Name: _____

MD/ENT Name: _____

PO#: _____

Bill to: _____

Ship to: _____

Shipping instructions: _____

03 Payment Method

I will pursue the Kanso 2 Sound Processor Upgrade Kit in partnership with my clinic and/or other funding means.

The self-pay discounted price does not apply.

I will self-pay the Kanso 2 Sound Processor Upgrade Kit and order directly through Cochlear.

I would like the CAD \$3,500 trade-in allowance to apply towards this upgrade order.¹

Allow one non-retired Nucleus Sound Processor trade-in for each Kanso 2 Upgrade Kit. Cannot be applied to the second sound processor option and this trade-in option is not available through Cochlear's Reimbursement & Insurance Services. Trade-in option must be declared at the time of the upgrade order and cannot be added after the upgrade has shipped.

I would like to keep in my existing non-retired processor.

04 Nucleus Kanso 2 Sound Processor

Select one for each item. Kit also includes: Home Charger, Magnet Removal Tool, Softwear Pads (pack of five), Headband (for pediatric recipients only), Socket Cover, Microphone Covers (pack of two), Cleaning Cloth, Screwdriver and Cleaning Brush, and User Guide.

Sound Processor Color:

Black Chocolate Brown Slate Grey
 Sandy Blonde Silver

Magnet Strength:

Profile™ Plus:

½(I) 1(I) 2(I) 3(I) 4(I) 5(I)

Standard:

½M 1M 2M 3M 4M 5M 6M

Not sure. Please use current magnet strength on file.

Safety Line:

Long (with Shirt Clip) Short (with Shirt Clip)
 Short Loop (with Hair Clips [pack of five])
 Select Hair Clip Color: White Brown Black

Kanso Halo Accessory:

An additional option to keep your Kanso 2 in place and is secured with hair clips on each side.

Kanso Halo Accessories (two, with Hair Clips [pack of five])
 Select Hair Clip Color: White Brown Black

Portable Charger Cable Length:

Short (12") Long (20")

05 Mix and Match Second Sound Processor²

Magnet strength will be the same strength selected in **SECTION 04** above.

I would like to choose a second sound processor for half off the self-pay price.

Kanso 2 Sound Processor

Includes: Sound Processor, Home Charger, Magnet, Magnet Removal Tool, Safety Line (same as selection made in Section 3), Softwear Pads (pack of five), Headband (for pediatric recipients only), Socket Cover, Microphone Covers (pack of two), Cleaning Cloth, Screwdriver and Cleaning Brush, and User Guide.

Sound Processor Color:

Black Chocolate Brown Slate Grey
 Sandy Blonde Silver

Kanso Halo Accessories (two, with Hair Clips [pack of five])

Select Hair Clip Color: White Brown Black

Nucleus 8 Sound Processor

Select one for each item below.

Kit includes: SoundProcessor, Power Extend Rechargeable Battery, Slimline Coil, Magnet and Earhook.

Sound Processor Color:

Black Grey Silver
 White Sand Brown

Slimline Integrated Coil Color:

Black Grey Silver
 White Sand Brown

Slimline Integrated Coil Cable Length:

(2" recommended)

2" (6cm) 3" (8cm) 4" (11cm) 10" (25cm)

Earhook Size (pack of three):

Small (Tamper Resistant) Small
 Medium Large

06 Make It Your Own

Select **two** options from below.

Kanso 2 Aqua+ Kit

A waterproof accessory that allows you to swim underwater, in salt water, and completely submerged.*

Includes: Aqua+, Activity Case, Headband (black) and Safety Line (Long).

Headband Size: XS (40–47cm) S (45–52cm) M (48–57cm) L (53–62cm)

Cochlear True Wireless Accessories

CR310 Remote Control

Recommended for recipients without an iPhone® or compatible Android™ smartphone.** Change programs and volume, and stream wirelessly.

Cochlear Wireless Phone Clip

Call hands-free and stream music from a Bluetooth® device.

Cochlear Wireless Mini Microphone 2+

Hear better over a distance or in noisy listening situations.

Cochlear Wireless TV Streamer

Stream high quality audio directly from your TV, stereo, computer, or other similar device.

07 Ready-to-Wear

Ready-to-Wear (Free of charge)

Have your Kanso 2 Sound Processor pre-loaded with your current programs by our Cochlear Hear Always Specialists. To participate, recipient must be age 5 or older. Cochlear recommends that you schedule a follow-up appointment with your audiologist to maximize your battery use and check your magnet strength. The Ready-to-Wear service may not be available for certain implant types. If your clinician decides that this option is not right for you, your Kanso 2 Sound Processor will be shipped unprogrammed to the address you provide.

08 Submit Form

Email completed order form to customer@cochlear.com

Policy:

Return and Exchange Policy: 90-Day Trial Period. 90 days to exchange unopened accessories.

Warranty: 3-years. No questions asked.

1. You may trade-in one (1) non-obsolete Nucleus Sound Processor when you choose to self-pay. This offer excludes upgrade orders placed through Reimbursement & Insurance Services. Limit 1 trade-in per upgraded device.
2. Second Sound Processor Option: Available within 90 days of upgrade order. No trade-in.

* The Kanso 2 Sound Processor with Aqua+ is dust and water resistant to level of IP68 of the International Standard IEC60529. The Kanso 2 Sound Processor with Aqua+ can be continuously submerged under water to a depth of up to 3 meters for up to 2 hours.

** For compatibility information visit www.cochlear.com/compatibility.

The Kanso 2 Sound Processor is compatible with the Nucleus Profile Plus, Nucleus Profile Series, CI500 Series, CI24RE (Freedom) Series, CI24R and CI24M implants. The Kanso 2 Sound Processor will not be compatible with the Nucleus® 22 (N22) Implant due to the Kanso 2 Sound Processor's size and coil type.

For more information on terms and conditions, please visit www.cochlear.com/us/upgrade.

Android and Google Play are registered trademarks of Google Inc.

©2023. Apple, the Apple logo, FaceTime, Made for iPad logo, Made for iPhone logo, Made for iPod logo, iPhone, iPad Pro, iPad Air, iPad mini, iPad and iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Cochlear is under license.

©Cochlear Limited 2023. All rights reserved. ACE, Advance OffStylet, AOS, Ardium, AutoNRT, Autosensitivity, Baha, Baha SoftWear, BCDrive, Beam, Bring Back the Beat, Button, Carina, Cochlear, 科利耳, コクレア, 코클리어, Cochlear SoftWear, Contour, 콘트واء, Contour Advance, Custom Sound, DermaLock, Freedom, Hear now. And always, Hugfit, Human Design, Hybrid, Invisible Hearing, Kanso, LowPro, MET, MP3000, myCochlear, mySmartSound, NRT, Nucleus, Osia, Outcome Focused Fitting, Off-Stylet, Piezo Power, Profile, Slimline, SmartSound, Softip, SoundArc, SoundBand, True Wireless, the elliptical logo, Vistafix, Whisper, WindShield and Xidium are either trademarks or registered trademarks of the Cochlear group of companies.

FUN3913 ISS7 AUG23

www.cochlear.com/us/upgrade • 800 587 6927 • customer@cochlear.com

Patient Name: _____