

Cochlear™ Osia® and Nucleus® System MRI Components Order Form

Before proceeding with an MRI, please confirm the specific implant type and refer to www.cochlear.us/mri for specific details and guidance about MRI Indications.

*** REQUIRED FIELDS**

01 Facility Information

Facility Name: _____	Bill to: _____
Facility Email: _____	Ship to: _____
Contact Name: _____	_____
PO#: _____	Shipping instructions: _____
Account #: _____	_____
	Date Required: _____
Additional comments (if needed): _____	

02 MRI Kits

Enter quantity required, up to five kits.

Includes: Round Magnetic Splints (two), Chinstraps (one Small, one Large), Bandage, Cochlear MRI Kit User Guide, for applicable product MRI Safety Checklist for Professionals and MRI Guidelines for Osia Implants.

By checking this box, you acknowledge that you are aware and in agreement that the Cochlear MRI Kit for Nucleus and Osia is intended to be applied by an Ear, Nose, and Throat (ENT) Physician.*

ACKNOWLEDGEMENT REQUIRED TO PROCESS ORDER

When applying the Cochlear Nucleus MRI Kit or Cochlear Osia MRI Kit, the following additional items may be needed:
Two hair ties, hairbrush and marker

___ **Cochlear Nucleus MRI Kit** P2025130

Note: Profile Plus Implant does not require the use of a Nucleus MRI Kit.

___ **Cochlear MRI Kit for Osia** P1890818

Note: OSI300 Implant does not require the use of an Osia MRI Kit.



03 Implant Magnet Removal Supplies

If indicated, enter quantity of item(s). Sterile products are single use only.

Nucleus System Temporary Implant Magnet Removal

CI600 Implant

- ___ Nucleus Non-Magnetic Cassette P782484
___ Nucleus Replacement Magnet Cassette P782485

CI500 Implant

- ___ CI500 Series Non-Magnetic Plug Z146624
___ CI500 Series Sterile Replacement Magnet Z179608

24RE Implant

- ___ CI24RE Non-Magnetic Plug Z50100
___ CI24RE Sterile Replacement Magnet Z50101

Retention item for use without Implant Magnet

- ___ Retainer Discs (50 pack) FUZ061

Osia System Temporary Implant Magnet Removal

OSI300 Implant

- ___ Cochlear Sterile Non-Magnetic Cassette P782484
___ Cochlear Sterile Replacement Magnet Cassette P1773917

OSI200 Implant (Prior Generation)

- ___ Cochlear Sterile Non-magnetic Plug P1620901
___ Cochlear Sterile Replacement Magnet P1620873

Long-term Use without Implant Magnet*

For maximum visibility, the Osia System may be used with the following items:

- ___ **OSI300 Implant:** Cochlear Sterile Non-Magnetic Cassette P782484 and Retainer Discs (50 pack) FUZ061
___ **OSI200 Implant:** Cochlear Sterile Non-magnetic Plug P1620901 and Retainer Discs (50 pack) FUZ061

04 Submit Form

Email or fax completed order form to procare@cochlear.com or 303 790 1157

Policy:

Returns: Unopened Bone Anchored Solution (BAS) and Nucleus system components can be returned within 30 days of the shipment date. Open BAS and Nucleus system components cannot be returned.

Exchanges: Unopened BAS System and Nucleus System components can be exchanged for the same generation technology within 90 days of the shipment date or fitting date for sound processors. Open Osia or Nucleus sound processors can be exchanged for the same generation technology within 90 days of fitting date. Exchanges are not accepted on any other open components.

* If a patient requires multiple MRI examinations at 3.0T or maximum visibility of the head is required, Osia may be used with a non-magnetic plug, and the Osia Sound Processor is retained by retainer disks provided by Cochlear. If serial imaging is no longer required, the non-magnetic plug may be replaced by a magnet for everyday usage. For details, please refer to our MRI Guidelines: www.cochlear.us/mri

©Cochlear Limited 2023. All rights reserved. ACE, Advance Off-Stylet, AOS, Ardiium, AutoNRT, Autosensitivity, Baha, Baha SoftWear, BCDrive, Beam, Bring Back the Beat, Button, Carina, Cochlear, 科利耳, コクレア, 코클리어, Cochlear SoftWear, Contour, 콘트우아, Contour Advance, Custom Sound, DermaLock, Freedom, Hear now. And always, Hugfit, Human Design, Hybrid, Invisible Hearing, Kanso, LowPro, MET, MP3000, myCochlear, mySmartSound, NRT, Nucleus, Osia, Outcome Focused Fitting, Off-Stylet, Piezo Power, Profile, Slimline, SmartSound, Softip, SoundArc, True Wireless, the elliptical logo, Vistafix, Whisper, WindShield and Xidium are either trademarks or registered trademarks of the Cochlear group of companies.

OSI132 ISS4 NOV23