

Cochlear[™] Kanso[®] 3 Nexa[®] Sound Processor (CP1175)

Features and components

- 1 Microphone ports
- 2 Indicator light
- 3 Magnet
- 4 Socket cover



Colour options



Black



Slate Grey



Silver



Sandy Blonde



Chocolate Brown

Technical specifications

Model number	CP1175
Wireless frequency	Proprietary 2.4 GHz low power bidirectional wireless link. Published commercial wireless protocol (Bluetooth® Low Energy)
Wireless communication range	Remote control: at least 2 m Smartphone control: at least 2 m Smartphone streaming: at least 7 m
Operating frequency	5.0 MHz
Input operating voltage	4.75 V to 5.35 V
Power consumption	15 mW–65 mW
Tap functions	Turn processor on and off, enable flight mode
Battery capacity	175 mAh (3.7V nominal)
Power options	Integrated rechargeable lithium ion battery (typical expected battery life 21 hours*)
Charge cycles	≥ 80% capacity after 2000 charge/discharge cycles at room temperature
Device dimensions	Length: 38 mm Width: 34 mm Depth: 12.5 mm Weight: 14.2 g [^]

Key features

- Smallest and lightest rechargeable off-the-ear sound processor
- SmartSound® iQ 2 with SCAN 2
- ForwardFocus (if enabled by a clinician) and SCAN 2 FF (automated ForwardFocus)[†]
- Dynamic power management
- Dual microphones
- Wireless inductive charging
- Integrated rechargeable lithium ion battery
- Available in five colour options
- IP rating for water and dust resistance:
 - IP68[~]
 - Swimmable with Aqua+
- Six magnet strengths
- LED status monitoring

Implant compatibility

- Compatible with Nucleus Nexa Series implants

Connectivity

- Ready for Bluetooth® LE Audio – firmware updates required[¥]
- Compatible with the Nucleus Smart App
- Direct streaming from compatible Apple and Android™ devices[#]
- Streaming from other Bluetooth-enabled devices via the Cochlear Wireless Phone Clip
- Compatible with Cochlear True Wireless™ devices
- Compatible with Cochlear Remote Control (CR310)

FM compatibility

- Via the Cochlear Wireless Mini Microphone 2+ with a three-pin Euro receiver

Accessories

- Retention accessories:
 - A range of safety lines
 - Kanso Headband (4 sizes)
- Cochlear SoftWear™ Pad
- Aqua+
- Kanso Halo accessory
- Portable Charger
- Home Charger

Programming

- Compatible with Custom Sound® fitting software
- Programming achieved using the programming pod, or wireless programming pod with the Kanso 2 programming adaptor cable

 Cochlear Ltd (ABN 96 002 618 073) 1 University Avenue, Macquarie University, NSW 2109, Australia T: +61 2 9428 6555

www.cochlear.com

* Battery life varies for every user, according to the age of the battery, the programs used each day, your implant type, the thickness of skin covering your implant, and the size and type of battery used. Streaming from compatible devices, True Wireless devices or FM may decrease sound processor battery life depending on how often and for how long streaming is engaged.

[^] Sound processor with lowest strength magnet

[†] ForwardFocus and SCAN 2 FF are clinician-enabled within the Custom Sound Pro fitting software and can be user-controlled with the Cochlear Nucleus Smart App or as an automated program.

[~] The Kanso 3 Sound Processor is dust and water resistant to level of IP68 of the International Standard IEC60529. With the Aqua+ accessory, the Kanso 3 Sound Processor is dust and water resistant to level of IP68 of the International Standard IEC60529 and can be continuously submerged under water to a depth of up to 3 metres for up to 2 hours. Refer to the relevant User Guide for more information.

[¥] As Bluetooth LE Audio compatible devices become available, a firmware update will be required for you to use certain features. Auracast™ broadcast audio capability is subject to third party adoption of the Auracast protocol.

[#] The Cochlear Nucleus Smart App is available on App Store and Google Play. The Kanso 3 Sound Processors are compatible with Apple and Android devices, for compatibility information visit www.cochlear.com/compatibility

Please seek advice from your health professional about treatments for hearing loss. Outcomes may vary, and your health professional will advise you about the factors which could affect your outcome. Always follow the directions for use. Not all products are available in all countries. Please contact your local Cochlear representative for product information. The product images shown are for illustrative purposes only and may not be an exact representation of the product. In Australia, Cochlear™ Nucleus® implant systems are intended for the treatment of moderately severe to profound hearing loss. For Cochlear™ implant systems: This product is not available for purchase by the general public. For information on funding and reimbursement please contact your healthcare professional.

Cochlear, 科利耳, コクレア, 코클리어, Hear now. And always, Nexa, Nucleus, Kanso, Advance Off-Stylet, AutoNRT, 콘트우아, Contour Advance, Custom Sound, Freedom, NRT, SmartSound, the elliptical logo, and marks bearing an ® or ™ symbol, are either trademarks or registered trademarks of the Cochlear group of companies (unless otherwise noted).

Android is a trademark of Google LLC. The Android robot is reproduced or modified from work created and shared by Google and used according to terms described in the Creative Commons 3.0 Attribution License. Google Play and the Google Play logo are trademarks of Google LLC.

Apple, the Apple logo, FaceTime, Made for iPad logo, Made for iPhone logo, Made for iPod logo, iPhone, iPad Pro, iPad Air, iPad mini, iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Cochlear Limited is under license.

© Cochlear Limited 2025. D2365655-V1 2025-07